# Xpert® FFPE Lysis Kit

REF

GXFFPE-LYSIS-CE-10

IVD

In Vitro Diagnostic  
Medical Device

©2020 Cepheid

301-5224 Rev. D, September 2020

Page 1

## English

*In Vitro* Diagnostic Medical Device.

## Description

The Xpert® FFPE Lysis Kit is designed to lyse Formalin-Fixed Paraffin Embedded (FFPE) tissue and preserve nucleic acids for subsequent GeneXpert analysis.

## Materials Provided

Xpert FFPE Lysis Kit (1 kit, sufficient to prepare 10 lysates)

Each kit contains:

- 1 bag containing 10 x 1.5 mL tubes
- 1 bag containing 10 x 5 mL vials
- 1 bottle containing 13 mL of bulk FFPE Lysis Reagent
- 1 tube containing 250 µL of bulk Proteinase K Reagent

## Materials and Equipment Required But Not Provided

- Hematoxylin and Eosin (H&E)
- Positively charged glass microscope slides
- Disposable Microtome Blades
- Microtome
- Water Bath (40 °C)

- Microscope slide drying rack
- 40-45 °C oven/incubator for drying slides (optional)
- Disposable razor blades (for macrodissection)
- Pipettes and filter pipette tips suitable to accurately transfer 5  $\mu$ L, 20  $\mu$ L, 260  $\mu$ L and 1.2 mL
- Heat Block suitable to hold 1.5 mL microcentrifuge tubes (and hold at 80 °C)
- Bench top vortex mixer
- Standard bench top microcentrifuge with fixed angle rotor that fits 1.5 mL microcentrifuge tubes
- $\geq$  95% Reagent Grade Ethanol
- Disposable gloves
- Adhesive labels or sample identification information

## Warning and Precautions

- All biological samples should be treated with standard universal precautions. Samples should only be handled by personnel trained in handling biohazardous materials.

- Handle all sample and kit reagents using appropriate techniques to prevent or minimize RNase and/or DNase contamination.

- Do not reuse macrodissection blades, pipette tips or tubes/vials to avoid cross contamination during sample handling.
- Spilled or leaking reagent tubes should be discarded and not used.
- Follow your institution's safety procedures for working with chemicals and handling biological samples.
- Incomplete removal (scraping) of tumor area from the slide for preparation of the FFPE lysate may result in insufficient material for the assay and therefore a higher than expected indeterminate/**INVALID** rate with Xpert Assays.

## Chemical Hazards

**Hazard Statement:** H319 Causes serious eye irritation

**Precautionary Statements:**

P264: Wash thoroughly after handling

P280: Wear protective gloves/protective clothing/eye protection/face protection

P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313: If eye irritation persists: Get medical advice attention.

## Kit Storage Requirements

Store the FFPE lysis kit at (2 °C to 28 °C). The expiration date for all lysis kit components is listed on the outer box labeling for all kit components.

## FFPE Tissue Requirements

1. Specimens must have been fixed in only 10% Neutral Buffered Formalin (NBF) for 6 to 72 hours to be tested with any Xpert assay that requires RNA or DNA from FFPE tissue.
2. Some FFPE tissue greater than 10 years old may be of insufficient quality for GeneXpert analysis.
3. The pathologist should select the FFPE tissue/tumor block with the greatest area of visible acceptance tissue/tumor for the assay to run.
4. The Xpert assay requires unstained slide mounted tissue for processing. If macrodissection is required, use an adjacent H&E stained slide from the FFPE tumor block as a guide to ensure the tumor area identified on the H&E stained slide is representative of the tumor area on the unstained slide.
5. Refer to the Xpert assay package insert for additional information.

## Slide/Scroll Preparation

Preparation of FFPE tissues using this lysis kit requires unstained slide mounted tissue section(s) or scroll(s) for processing and an adjacent H&E stained slide from the same tissue block.

1. Using a microtome, cut a 4-5  $\mu\text{m}$  thick section/scroll for H&E staining.
2. Using a microtome, cut a sufficient number of 4-5  $\mu\text{m}$  thick sections/scrolls adjacent to the section designated for H&E, as determined by pathologist, for use in the FFPE Lysis kit.  
For specimens that require multiple slides/scrolls to meet minimum Xpert assay requirements, all slides/scrolls must be processed together.
3. To process tissue sections follow the sections below.
  - a. Float the section(s) in a water bath at 40 °C.
  - b. Mount the section(s), one per slide, onto positively charged glass microscope slides.
  - c. Allow slide(s) to air dry overnight or bake at 40-45 °C overnight (optional).
  - d. Slide section(s) should be stored at 2 °C to 8 °C and processed within 1 month after preparation. Freshly cut sections or unstained slides will yield highest quality and most reliable results as mRNA degrades more rapidly on slide sections exposed to air than in FFPE blocks.
4. To process tissue scroll(s) follow the steps below.
  - a. Place scroll(s) in a provided 1.5 mL lysis tube. Label the tube for each sample to be processed.
  - b. Scroll(s) in a 1.5 mL tube should be stored at 2 °C to 8 °C and processed within 2 weeks after preparation.

## Tissue Removal from Slide

1. Label a 1.5 mL lysis tube (provided) for each sample to be processed.
2. If macrodissection is not required:
  - a. Using a new razor blade or scalpel for each tissue sample to be processed, completely remove (scrape) the invasive tumor tissue section from the slide and transfer to the labeled 1.5 mL lysis tube.
  - b. Section(s) in a 1.5 mL tube should be stored at 2°C to 8°C and processed within 2 weeks.
3. If macrodissection is required:
  - a. Examine the H&E stained slide (pathologist).  
Identify (and outline) the tumor area for the assay. Refer to Xpert assay package insert for required number of slides or minimum tumor cellularity requirement.
  - b. Prepare for macrodissection (pathologist).  
Outline the tumor area to be used for the assay on the backside of the unstained slide(s) by aligning it with the corresponding H&E stained slide and transposing the outlined area.
  - c. Perform macrodissection (pathologist or technician).  
Using a new razor blade or scalpel for each tissue sample to be processed, completely remove (scrape) the outlined invasive tumor tissue from the slide (see Figure 1) and transfer to the labeled 1.5 mL lysis tube.

**Figure 1: Examples of Proper (Recommended) and Improper (Not Recommended) Tissue Removal from Slide**

- Macrodissected section(s) in a 1.5 mL tube should be stored at 2°C to 8°C and processed within 2 weeks.

## FFPE Tissue Processing

1. Preheat heat block to 80 °C.
2. Add 1.2 mL of FFPE lysis reagent to the 1.5 mL tube containing FFPE section/scroll.
3. Add 20 µL of Proteinase K (PK) to the same 1.5 mL lysis tube.
4. Close lid.
5. Vortex the sample continuously at a maximum setting for 5 seconds.
6. Briefly microcentrifuge the tube to remove liquid from the lid.
7. Using a heat block set at 80 °C, incubate 1.5 mL lysis tube containing sample and lysis reagent for 30 minutes.
8. Vortex the sample continuously at maximum setting for 5 seconds.
9. Briefly microcentrifuge the tube to remove liquid from the lid.
10. Using a pipette, transfer the entire contents (~1.2 mL) to a provided 5 mL sample vial.
11. Label the vial for each sample to be processed.
12. Add 1.2 mL of ≥ 95% ethanol to the same 5 mL sample vial.
13. Secure the cap and vortex the sample continuously at maximum setting for 15 seconds.


## FFPE Tissue Processing – Concentrated Lysate

Follow this protocol if a more concentrated lysate is desired due to insufficient sample (This should be considered for  $\leq 6 \times 1 \text{ mm}^2$  area of tumor on scraped slide). Refer to Xpert assay package insert for additional information.

1. Preheat heat block to  $80^\circ\text{C}$ .
2. Add  $260 \mu\text{L}$  of FFPE lysis reagent to the 1.5 mL lysis tube containing the FFPE section.
3. Add  $5 \mu\text{L}$  of Proteinase K (PK) to the same 1.5 mL lysis tube.
4. Close lid.
5. Vortex the sample continuously at maximum setting for 5 seconds.
6. Briefly microcentrifuge the tube to remove liquid from the lid.
7. Using a heat block set at  $80^\circ\text{C}$ , incubate the 1.5 mL lysis tube containing sample and lysis reagent for 30 minutes.
8. Vortex the sample continuously at maximum setting for 5 seconds.
9. Briefly microcentrifuge the tube to remove liquid from the lid.
10. Add  $260 \mu\text{L}$  of  $\geq 95\%$  ethanol to the same 1.5 mL lysis tube.
11. Close lid.
12. Vortex the lysis tube continuously at maximum setting for 15 seconds.
13. Briefly microcentrifuge the tube to remove liquid from the lid.

## Sample Storage and Transport

 Prepared lysate, with ethanol, should be transported to the laboratory at  $2^\circ\text{C}$  to  $8^\circ\text{C}$ , if testing is to be performed within 1 week. If assay is to be performed at a later time, the FFPE lysate is stable and may be stored up to 4 weeks at  $\leq -20^\circ\text{C}$  before testing.

## Table of Symbols

Symbol	Meaning
	Catalog number
	CE Marking - European Conformity
	Contains sufficient for <n> tests
	Manufacturer
	Country of manufacture
	Temperature limitation
	Batch code
	Consult instructions for use
	Biological Risks
	<i>In Vitro</i> Diagnostic Medical Device
	Authorized Representative in the European Community
	Do not reuse
	Caution
	Warning
	Expiration Date

## Technical Assistance

<b>United States</b>	<b>France</b>
Telephone: + 1 888 838 3222	Telephone: + 33 563 825 319
Email: techsupport@cepheid.com	Email: support@cepheideurope.com

Contact information for all Cepheid Technical Support offices is available on our website: [www.cepheid.com/en/CustomerSupport](http://www.cepheid.com/en/CustomerSupport).

Manufactured for Cepheid  
Cepheid  
904 Caribbean Drive  
Sunnyvale, CA 94089 USA  
Phone: +1.408.541.4191  
Fax: +1.408.541.4192

Cepheid Europe SAS  
Vira Solelh  
81470 Maurens-Scopont  
France  
Phone: +33.563.825.300  
Fax: +33.563.825.301

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