## English

For *in vitro* diagnostic use only.

### Intended Use

The Xpert® Swab Specimen Collection Kit is designed to collect, preserve, and transport endocervical, vaginal, pharyngeal, and rectal swab specimens. This transport system is used for testing with Xpert tests.

### Summary and Explanation of the Procedure

Endocervical (ES), vaginal (VS), pharyngeal (PS), or rectal (RS) swab specimens are collected from patients using flocked swabs included in the specimen collection kit. Swabs are broken off into the transport reagent tubes to elute organisms and stabilize DNA. Swab specimens are then transported to the laboratory for testing on the GeneXpert Instrument Systems.

### Materials Provided

Xpert Swab Specimen Collection Kits: (quantity: 50)

Each kit contains a sealed package containing:

- 1 individually wrapped sterile large cleaning swab
- 1 package containing:
  - 1 tube containing Xpert Swab Transport Reagent (pink cap)
  - 1 individually wrapped sterile flocked collection swab
- Instructions for collection of patient-collected vaginal swabs (quantity: 50)
- Instructions for collection of endocervical, pharyngeal, and rectal swabs (quantity: 1 of each per kit)

### Materials Required but Not Provided

- Disposable gloves
- Adhesive labels or permanent markers for recording specimen information

### Swab Specimens Appropriate for Use

The table below identifies swab specimen types and targets that are acceptable for use with the Xpert CT/NG and the Xpert TV tests.

Xpert Test	Target	ES	VS	PS	RS
Xpert CT/NG	<i>Chlamydia trachomatis</i> and <i>Neisseria gonorrhoeae</i>	Yes	Yes	Yes	Yes
Xpert TV	<i>Trichomonas vaginalis</i>	Yes	Yes	No	No

ES = Endocervical Swab, VS = Vaginal Swab,  
PS = Pharyngeal Swab, RS = Rectal Swab

### Warnings and Precautions

- For prescription use only.
- Treat all biological specimens, including used cartridges, as if capable of transmitting infectious agents. Because it is often impossible to know which might be infectious, all biological specimens should be treated with standard precautions. Guidelines for specimen handling are available from the U.S. Center for Disease Control and Prevention and the Clinical and Laboratory Standards Institute.<sup>1,2</sup>
- Avoid cross contamination during sample handling steps.
- Do not use kit if opened or damaged.
- Spilled or leaking transport reagent tubes should be discarded and not be used.

- Use this collection kit only with the Xpert CT/NG or the Xpert TV tests.
- Follow your institution's safety procedures for working with chemicals and handling biological samples.
- With vaginal/endocervical specimens tested with the Xpert CT/NG test, interference may be observed in the presence of blood (>1% v/v) or mucin (>0.8% w/v).
- With vaginal/endocervical specimens tested with the Xpert TV test, interference may be observed in the presence of blood (>60% v/v).

### Kit Storage Requirements

Store collection kit at 2 °C to 30 °C.

**IMPORTANT:** Do not use after the expiration date which is clearly printed on the outer box label.

### Product Deterioration

The Xpert Swab Collection Kit should not be used if (1) there is evidence of damage or contamination to the product, (2) there is evidence of leakage, (3) the expiration date has passed on any of the components, or (4) the swab pouch is open.

### Limitations

- Specimens should be handled aseptically.
- The Xpert Swab Collection Kit is intended to be used with the transport reagent tube and swab provided in the kit. The use of transport reagent tubes or swabs from any other source could affect the performance of the product.
- Use this collection kit only with the Xpert CT/NG or the Xpert TV tests.

### Specimen Collection Procedure Patient-collected Vaginal Swab Specimen Collection

**Caution:** Do NOT expose swab to Xpert Swab Transport Reagent prior to collection.

**Caution:** Wash your hands before starting.

**Note:** If you have any questions about this procedure, please ask your doctor, nurse, or care-provider.

1. In the privacy of the examination room or restroom, undress from the waist down. Comfortably position yourself to maintain balance during the collection procedure.
2. Open the outer peelpack (which contains the two-package kit), and identify the larger cleaning swab and discard it.
3. Open the package that contains the pink-capped Xpert Swab Transport Reagent tube and individually wrapped collection swab. Set the tube aside before beginning to collect sample.
4. Open the collection swab wrapper by peeling open the top of the wrapper.
5. Remove the swab, taking care not to touch the tip or lay it down. If the soft tip is touched, the swab is laid down, or the swab is dropped, request a new collection kit.
6. Hold the swab in your hand, placing your thumb and forefinger in the middle of the swab shaft across the scoreline.
7. Carefully insert the swab into your vagina about 5 cm (two inches) inside the opening of the vagina and gently rotate the swab for 10 to 30 seconds. Ensure the swab touches the walls of the vagina so that moisture is absorbed by the swab.
8. Withdraw the swab carefully and continue to hold it in your hand.
9. While holding the swab in the same hand, unscrew the pink cap from the Xpert Swab Transport Reagent tube.

10. Do not spill the contents of the tube. If the contents of the tube are spilled, request a new collection kit.

**WARNING:** If the contents of the tube are spilled on your skin, wash the affected area with soap and water. If the contents of the tube are splashed in your eyes, immediately flush your eyes with water. Notify your doctor, nurse or care-provider if irritation develops. If the contents of the tube are spilled, your test result may be invalidated. Do not take internally.

11. Immediately place the specimen collection swab into the transport reagent tube.
12. Identify the scoreline on the collection swab shaft. Carefully break the swab shaft against the side of the tube at the scoreline and discard the top portion of the swab shaft. If needed, gently rotate the swab shaft to complete the breakage. Use care to avoid splashing contents on the skin. Wash with soap and water if exposed.
13. Re-cap the swab transport reagent tube and tighten the cap securely.
14. Return the tube as instructed by your doctor, nurse, or care-provider to complete the following steps.

### Steps For Caregiver To Perform

1. Invert or gently shake the tube 3-4 times to elute material from the swab. Avoid foaming.
2. Label the transport tube with the sample identification information, including date of collection, as required.
3. Transport and store the swab samples in the Xpert Swab Transport Reagent tube according to the conditions as indicated in Table 1 for the Xpert CT/NG test or for the Xpert TV test (as appropriate).

### Clinician-collected Endocervical Swab Specimen Collection

**Caution:** Do NOT expose swab to Xpert Swab Transport Reagent prior to collection.

1. Open the Xpert Swab Specimen collection kit.
2. Before collecting the endocervical specimen with the Xpert Swab Specimen Collection Kit, remove excess mucus from the cervical os and surrounding mucosa using the large individually wrapped cleaning swab. Partially peel open the cleaning swab wrapper and remove the swab. Clean the cervical os and surrounding mucosa and then discard the swab.

**Note:** If collecting multiple specimens, excess mucus need only be removed once.

3. Open the package that contains the pink-capped Xpert Swab Transport Reagent tube and the individually wrapped collection swab. Set the tube aside before proceeding.
4. Open the collection swab wrapper by peeling open the top of the wrapper.
5. Remove the swab, taking care not to touch the tip or lay it down. If the soft tip is touched, the swab is laid down, or the swab is dropped, use a new Xpert Swab Specimen Collection Kit.
6. Insert the collection swab into the endocervical canal.
7. Gently rotate the swab clockwise for 10 to 30 seconds in the endocervical canal to ensure adequate sampling.
8. Withdraw the swab carefully.
9. While holding the swab in the same hand, unscrew the cap from the Xpert Swab Transport Reagent tube.

10. Do not spill the contents of the tube. If the contents of the tube are spilled, use a new collection kit.

**WARNING:** If the contents of the tube are spilled on your skin, wash the affected area with soap and water. If the contents of the tube are splashed in your eyes, immediately flush your eyes with water.

11. Immediately place the swab into the transport reagent tube.
12. Identify the scoreline on the collection swab shaft. Carefully break the swab shaft against the side of the tube at the scoreline and discard the top portion of the swab shaft; use care to avoid splashing the contents. Wash with soap and water if exposed.
13. Re-cap the swab transport reagent tube and tighten the cap securely.
14. Invert or gently shake the tube 3-4 times to elute material from the swab. Avoid foaming.
15. Label the transport tube with the sample identification information, including date of collection, as required.
16. Transport and store the swab samples in the Xpert Swab Transport Reagent tube according to the conditions indicated in Table 1 for the Xpert CT/NG test or for the Xpert TV test (as appropriate).

### Clinician-collected Pharyngeal Swab Specimen Collection

**Caution:** Do NOT expose swab to Xpert Swab Transport Reagent prior to collection.

1. Open the Xpert Swab Specimen collection kit.
2. Open the outer peelpack (which contains the two-package kit), and identify the larger cleaning swab and discard it.
3. Open the collection swab wrapper by partially peeling open the top of the wrapper.
4. Remove the swab, taking care not to touch the tip or lay it down. If the soft tip is touched, the swab is laid down, or the swab is dropped, use a new Xpert Swab Specimen Collection Kit.
5. Instruct the patient to open mouth widely. Position the tongue toward the bottom of the mouth. Swab areas of the pharynx (tonsil, posterior wall, uvula, posterior wall).
6. While holding the swab in the same hand, unscrew the cap from the Xpert Swab Transport Reagent tube.
7. Do not spill the contents of the tube. If the contents of the tube are spilled, use a new collection kit.

**WARNING:** If the contents of the tube are spilled on your skin, wash the affected area with soap and water. If the contents of the tube are splashed in your eyes, immediately flush your eyes with water.

8. Immediately place the swab into the transport reagent tube.
9. Identify the scoreline on the collection swab shaft. Carefully break the swab shaft against the side of the tube at the scoreline and discard the top portion of the swab shaft; use care to avoid splashing the contents. Wash with soap and water if exposed.
10. Re-cap the swab transport reagent tube and tighten the cap securely.
11. Invert or gently shake the tube 3-4 times to elute material from the swab. Avoid foaming.
12. Label the transport tube with the sample identification information, including date of collection, as required.
13. Transport and store the swab samples in the Xpert Swab Transport Reagent tube according to the conditions indicated in Table 1 for the Xpert CT/NG test.

## Clinician-collected Rectal Swab Specimen Collection

**Caution:** Do NOT expose swab to Xpert Swab Transport Reagent prior to collection.

1. Open the Xpert Swab Specimen collection kit.
2. Open the outer peelpack (which contains the two-package kit), and identify the larger cleaning swab and discard it.
3. Open the collection swab wrapper by partially peeling open the top of the wrapper.
4. Remove the swab, taking care not to touch the tip or lay it down. If the soft tip is touched, the swab is laid down, or the swab is dropped, use a new Xpert Swab Specimen Collection Kit.

**WARNING:** Do not wet or use lubricant on swab.

5. Carefully insert the swab approximately 1 cm beyond the anal sphincter (so that the fiber tips are no longer visible), and rotate gently for 10 seconds. See Figure 1 for examples of acceptable swab specimens to be used with Xpert CT/NG. See Figure 2 for examples of highly soiled swab specimens that should not be used with Xpert CT/NG.


Figure 1. Examples of Acceptable Rectal Swabs


Figure 2. Examples of Unacceptable Rectal Swabs

6. While holding the swab in the same hand, unscrew the cap from the Xpert Swab Transport Reagent tube.
7. Do not spill the contents of the tube. If the contents of the tube are spilled, use a new collection kit.

**WARNING:** If the contents of the tube are spilled on your skin, wash the affected area with soap and water. If the contents of the tube are splashed in your eyes, immediately flush your eyes with water.

8. Immediately place the swab into the transport reagent tube.
9. Identify the scoreline on the collection swab shaft. Carefully break the swab shaft against the side of the tube at the scoreline and discard the top portion of the swab shaft; use care to avoid splashing the contents. Wash with soap and water if exposed.
10. Re-cap the swab transport reagent tube and tighten the cap securely.
11. Invert or gently shake the tube 3-4 times to elute material from the swab. Avoid foaming.
12. Label the transport tube with the sample identification information, including date of collection, as required.
13. Transport and store the swab samples in the Xpert Swab Transport Reagent tube according to the conditions indicated in Table 1 for the Xpert CT/NG test.

Table 1. Specimen Transport and Storage

Specimen	Xpert Test	Transport and Storage Temperature (°C)	Storage Time
Endocervical Swab, Vaginal Swab, Pharyngeal Swab, and Rectal Swab in Xpert Swab Transport Reagent	Xpert CT/NG	2-30 °C	60 Days
Endocervical Swab and Vaginal Swab in Xpert Swab Transport Reagent	Xpert TV	2-30 °C	60 Days

## References

1. Centers for Disease Control and Prevention. Biosafety in microbiological and biomedical laboratories (refer to latest edition). [www.cdc.gov/biosafety/publications](http://www.cdc.gov/biosafety/publications).
2. Clinical and Laboratory Standards Institute. Protection of laboratory workers from occupationally acquired infections; Approved Guideline. Document M29 (refer to latest edition).

## Table of Symbols

Symbol	Meaning
	Catalog number
	CE marking – European Conformity
	<i>In vitro</i> diagnostic medical device
	Contains sufficient for <n> tests
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Authorized representative in the European Community
	Do not reuse
	Manufacturer
	Temperature limitation
	Batch code
	Caution
	For prescription use only
	Consult instructions for use
	Expiration date
	Biological risks
	Do not use if package is damaged
	Country of manufacture

## Technical Support Contact

### Contact Information

United States

Telephone: + 1 888 838 3222

Email: techsupport@cepheid.com

France

Telephone: + 33 563 825 319

Email: support@cepheideurope.com

Contact information for all Cepheid Technical Support offices is available on our website:

[www.cepheid.com/en/CustomerSupport](http://www.cepheid.com/en/CustomerSupport).

Cepheid  
904 Caribbean Drive  
Sunnyvale, CA 94089  
USA  
Phone: +1.408.541.4191  
Fax: +1.408.541.4192

Cepheid Europe SAS  
Vira Solelh  
81470 Maurens-Scopont  
France  
Phone: +33 563 825 300  
Fax: +33 563 825 301

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### Bibliographie

- Centers for Disease Control and Prevention. Biosafety in microbiological and biomedical laboratories (consulter l'édition la plus récente). [www.cdc.gov/biosafety/publications](http://www.cdc.gov/biosafety/publications).
- Clinical and Laboratory Standards Institute. Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline. Document M29 (consulter l'édition la plus récente).

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1. Centers for Disease Control and Prevention. Biosafety in microbiological and biomedical laboratories (consultar la última edición). [www.cdc.gov/biosafety/publications](http://www.cdc.gov/biosafety/publications).
2. Clinical and Laboratory Standards Institute. Protection of laboratory workers from occupationally acquired infections; Approved Guideline. Documento M29 (consultar la última edición).

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### Referências

1. Centers for Disease Control and Prevention. Biosafety in Microbiological and Biomedical laboratories (consultar a edição mais recente). [www.cdc.gov/biosafety/publications](http://www.cdc.gov/biosafety/publications).
2. Clinical and Laboratory Standards Institute. Protection of Laboratory Workers from Occupationally Acquired Infections; Guideline aprovada. Documento M29 (consulte a edição mais recente).

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- Centers for Disease Control and Prevention. Biosafety in microbiological and biomedical laboratories (см. последнюю редакцию). [www.cdc.gov/biosafety/publications](http://www.cdc.gov/biosafety/publications).
- Clinical and Laboratory Standards Institute. Protection of laboratory workers from occupationally acquired infections; Approved Guideline. документ M29 (см. последнюю редакцию).

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- Centers for Disease Control and Prevention. Biosafety in microbiological and biomedical laboratories (див. останнє видання). [www.cdc.gov/biosafety/publications](http://www.cdc.gov/biosafety/publications).
- Clinical and Laboratory Standards Institute. Protection of laboratory workers from occupationally acquired infections; Approved Guideline. Document M29 (див. останнє видання).

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