English

For In Vitro Diagnostic Use Only.

Description

The GeneXpert® Urine Transport Reagent Kit is designed to collect, preserve, and transport urine samples to the laboratory for analysis.

Materials Provided

GeneXpert Urine Transport Reagent Kit (quantity: 30)

Each kit contains:

- 1 disposable transfer pipette
- 1 tube containing 4.5 mL of Urine Transport Reagent (contains guanidinium chloride)

Materials Required but Not Provided

- Clean, plastic, preservative-free, urine specimen collection cups
- Disposable gloves
- Adhesive labels for sample identification information

Warnings and Precautions

- All biological specimens should be treated with standard precautions. Samples should only be handled by personnel trained in handling infectious materials.
- Avoid cross contamination during sample handling steps.
- Spilled or leaking transport reagent tubes should be discarded and not used.
- Follow your institution's safety procedures for working with chemicals and handling biological samples.
- Consult your institution's environmental waste personnel on proper disposal of unused reagents.
- Check state and local regulations as they may differ from federal disposal regulations. This material may exhibit characteristics of hazardous waste requiring specific disposal requirements. Institutions should check their country hazardous waste disposal requirements.
- Urine Transport Reagent contains guanidinium chloride which is toxic if swallowed (H302), irritating to eyes and skin (H316, H320), and is harmful to aquatic life (H402). See Regulation EC 1272/2008 (EU CLP).

Kit Storage Requirements

- Store the collection kit at room temperature (15–28 °C).
- Do not use reagents or cartridges that have passed the expiration date.

Sample Collection and Handling

1. Ensure that you have at least 4.5 mL of urine and that it is not from the first void of the day.
2. Invert the urine collection cup three times to mix.
3. Open one transport kit.
4. Remove the cap from the tube.
5. Remove the cap from the urine collection cup.
6. Find the black “fill to” line (9 mL mark) on the tube.
7. Use the transfer pipette to add urine to the “fill to” line.
8. Replace the tube cap and tighten.
9. Invert the tube three times to mix.

Sample Storage and Transport

Urine samples should be transferred to the GeneXpert Urine Transport Reagent tubes within 1 hour of sample collection. Urine samples stored in GeneXpert Urine Transport Reagent tubes should be transported to the laboratory at 2–28 °C. Urine samples in the Urine Transport Reagent tubes are stable up to 7 days at 2–28 °C before testing.

Table of Symbols

| Symbol | Meaning |
|---|---|
|  | Catalog number |
|  | <i>In vitro</i> diagnostic medical device |
|  | CE marking – European Conformity |
|  | Authorized Representative in the European Community |
|  | Contains sufficient for <n> tests |
|  | Do not reuse |
|  | Manufacturer |
|  | Temperature limitation |
|  | Batch code |
|  | Expiration date |
|  | Consult instructions for use |
|  | Date of manufacture |
|  | Warning |
|  | Biological risks |

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Cepheid Europe SAS