# Xpert<sup>®</sup> MTB/XDR

**REF** GXMTB/XDR-10

Petunjuk Penggunaan

**IVD** CE

## **Pernyataan Merek Dagang, Paten, dan Hak Cipta**

Cepheid<sup>®</sup>, the Cepheid logo, GeneXpert<sup>®</sup>, and Xpert<sup>®</sup> are trademarks of Cepheid, registered in the U.S. and other countries. All other trademarks are the property of their respective owners.

THE PURCHASE OF THIS PRODUCT CONVEYS TO THE BUYER THE NON-TRANSFERABLE RIGHT TO USE IT IN ACCORDANCE WITH THESE INSTRUCTIONS FOR USE. NO OTHER RIGHTS ARE CONVEYED EXPRESSLY, BY IMPLICATION OR BY ESTOPPEL. FURTHERMORE, NO RIGHTS FOR RESALE ARE CONFERRED WITH THE PURCHASE OF THIS PRODUCT.

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| Test Result | Analyte Result  | Detail  | Melt Peaks | Errors | History | Support |
|-------------|---|---------|------------|--------|---------|---------|
| Assay Name  | MTB-XDR   | Version | 3          |        |         |         |
| Test Result | <b>MTB DETECTED;</b><br><b>INH Resistance DETECTED;</b><br><b>FLQ Resistance DETECTED;</b><br><b>AMK Resistance DETECTED;</b><br><b>KAN Resistance DETECTED;</b><br><b>CAP Resistance DETECTED;</b><br><b>ETH Resistance DETECTED</b> |         |            |        |         |         |

| Test Result     | Analyte Result        | Detail           | Melt Peaks | Errors | History | Support |
|-----------------|-----------------------|------------------|------------|--------|---------|---------|
| Analyte Name    | Melt Peak Temperature | Melt Peak Height |            |        |         |         |
| inhA-melt       |                       |                  |            |        |         |         |
| katG-melt       |                       |                  |            |        |         |         |
| fabG1-melt      |                       |                  |            |        |         |         |
| ahpC-melt       |                       |                  |            |        |         |         |
| gyrA1-melt      | 76.1                  | 90.0             |            |        |         |         |
| gyrA2-melt      | 69.6                  | 39.7             |            |        |         |         |
| gyrA3-melt      |                       |                  |            |        |         |         |
| gyrB2-melt      |                       |                  |            |        |         |         |
| rrs-melt        |                       |                  |            |        |         |         |
| eis-melt        |                       |                  |            |        |         |         |
| inhA-mut melt   | 70.9                  | 259.6            |            |        |         |         |
| katG-mut melt   | 68.4                  | 214.0            |            |        |         |         |
| fabG1-mut melt  | 75.9                  | 181.1            |            |        |         |         |
| ahpC-mut melt   | 66.2                  | 68.2             |            |        |         |         |
| gyrA1-mutA melt |                       |                  |            |        |         |         |
| gyrA1-mutB melt |                       |                  |            |        |         |         |
| gyrA1-mutC melt |                       |                  |            |        |         |         |
| gyrA2-mutA melt |                       |                  |            |        |         |         |
| gyrA2-mutB melt |                       |                  |            |        |         |         |
| gyrA3-mutA melt |                       |                  |            |        |         |         |
| gyrA3-mutB melt | 76.0                  | 125.0            |            |        |         |         |
| gyrA3-mutC melt |                       |                  |            |        |         |         |
| gyrB2-mut melt  | 66.0                  | 103.2            |            |        |         |         |
| rrs-mut melt    | 71.0                  | 125.7            |            |        |         |         |
| eis-mutA melt   | 71.4                  | 163.9            |            |        |         |         |
| eis-mutB melt   |                       |                  |            |        |         |         |

**Gambar 9. MTB TERDETEKSI (MTB DETECTED); Resistansi INH, FLQ, AMK, KAN, CAP, dan ETH TERDETEKSI (DETECTED)**

| Test Result                  | Analyte Result   | Detail           | Melt Peaks | Errors | History | Support |
|------------------------------|--|------------------|------------|--------|---------|---------|
| Assay Name MTB-XDR Version 3 |  |                  |            |        |         |         |
| Test Result                  | <b>MTB DETECTED;</b><br><b>INH Resistance DETECTED;</b><br>FLQ Resistance NOT DETECTED;<br>AMK Resistance NOT DETECTED;<br>KAN Resistance NOT DETECTED;<br>CAP Resistance NOT DETECTED;<br>ETH Resistance NOT DETECTED |                  |            |        |         |         |
| Test Result                  | Analyte Result   | Detail           | Melt Peaks | Errors | History | Support |
| Analyte Name                 | Melt Peak Temperature  | Melt Peak Height |            |        |         |         |
| inhA-melt                    | 76.6   | 284.9            |            |        |         |         |
| katG-melt                    | 74.0   | 105.2            |            |        |         |         |
| fabG1-melt                   |  |                  |            |        |         |         |
| ahpC-melt                    | 69.0   | 35.4             |            |        |         |         |
| gyrA1-melt                   | 76.6   | 65.2             |            |        |         |         |
| gyrA2-melt                   | 70.4   | 64.9             |            |        |         |         |
| gyrA3-melt                   | 71.4   | 92.2             |            |        |         |         |
| gyrB2-melt                   | 69.7   | 84.7             |            |        |         |         |
| rrs-melt                     | 75.3   | 146.8            |            |        |         |         |
| eis-melt                     | 68.7   | 124.2            |            |        |         |         |
| inhA-mut melt                |  |                  |            |        |         |         |
| katG-mut melt                |  |                  |            |        |         |         |
| fabG1-mut melt               | 75.9   | 178.0            |            |        |         |         |
| ahpC-mut melt                |  |                  |            |        |         |         |
| gyrA1-mutA melt              |  |                  |            |        |         |         |
| gyrA1-mutB melt              |  |                  |            |        |         |         |
| gyrA1-mutC melt              |  |                  |            |        |         |         |
| gyrA2-mutA melt              |  |                  |            |        |         |         |
| gyrA2-mutB melt              |  |                  |            |        |         |         |
| gyrA3-mutA melt              |  |                  |            |        |         |         |
| gyrA3-mutB melt              |  |                  |            |        |         |         |
| gyrA3-mutC melt              |  |                  |            |        |         |         |
| gyrB2-mut melt               |  |                  |            |        |         |         |
| rrs-mut melt                 |  |                  |            |        |         |         |
| eis-mutA melt                |  |                  |            |        |         |         |
| eis-mutB melt                |  |                  |            |        |         |         |

Gambar 10. MTB TERDETEKSI (MTB DETECTED); Resistansi INH TERDETEKSI (DETECTED)

Test Result Analyte Result Detail Melt Peaks Errors History Support

Assay Name MTB-XDR Version 4

Test Result

MTB DETECTED;

INH Resistance DETECTED;

FLQ Resistance NOT DETECTED;

AMK Resistance INDETERMINATE;

KAN Resistance DETECTED;

CAP Resistance INDETERMINATE;

ETH Resistance DETECTED

Test Result Analyte Result Detail Melt Peaks Errors History Support

| Analyte Name    | Melt Peak Temperature | Melt Peak Height |
|-----------------|-----------------------|------------------|
| inhA-melt       |                       |                  |
| katG-melt       |                       |                  |
| fabG1-melt      | 71.5                  | 254.6            |
| ahpC-melt       | 68.7                  | 49.4             |
| gyrA1-melt      | 76.3                  | 62.9             |
| gyrA2-melt      | 70.2                  | 59.8             |
| gyrA3-melt      | 71.5                  | 56.5             |
| gyrB2-melt      | 69.4                  | 74.8             |
| rrs-melt        |                       |                  |
| eis-melt        |                       |                  |
| inhA-mut melt   | 70.9                  | 277.7            |
| katG-mut melt   | 68.2                  | 157.7            |
| fabG1-mut melt  |                       |                  |
| ahpC-mut melt   |                       |                  |
| gyrA1-mutA melt |                       |                  |
| gyrA1-mutB melt |                       |                  |
| gyrA1-mutC melt |                       |                  |
| gyrA2-mutA melt |                       |                  |
| gyrA2-mutB melt |                       |                  |
| gyrA3-mutA melt |                       |                  |
| gyrA3-mutB melt |                       |                  |
| gyrA3-mutC melt |                       |                  |
| gyrB2-mut melt  |                       |                  |
| rrs-mut melt    |                       |                  |
| eis-mutA melt   |                       |                  |
| eis-mutB melt   | 62.6                  | 46.5             |

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| Test Result   | Analyte Result        | Detail           | Melt Peaks | Errors | History | Support |
|---|-----------------------|------------------|------------|--------|---------|---------|
| Assay Name MTB-XDR Version 3  |                       |                  |            |        |         |         |
| <b>Test Result</b><br><b>MTB DETECTED;</b><br><b>INH Resistance DETECTED;</b><br><b>Low FLQ Resistance DETECTED;</b><br>AMK Resistance NOT DETECTED;<br>KAN Resistance NOT DETECTED;<br>CAP Resistance NOT DETECTED;<br>ETH Resistance NOT DETECTED |                       |                  |            |        |         |         |
| Test Result   | Analyte Result        | Detail           | Melt Peaks | Errors | History | Support |
| Analyte Name  | Melt Peak Temperature | Melt Peak Height |            |        |         |         |
| inhA-melt   | 76.5                  | 313.1            |            |        |         |         |
| katG-melt   |                       |                  |            |        |         |         |
| fabG1-melt  | 71.7                  | 211.5            |            |        |         |         |
| ahpC-melt   | 69.0                  | 47.2             |            |        |         |         |
| gyrA1-melt  |                       |                  |            |        |         |         |
| gyrA2-melt  |                       |                  |            |        |         |         |
| gyrA3-melt  |                       |                  |            |        |         |         |
| gyrB2-melt  | 69.6                  | 81.1             |            |        |         |         |
| rrs-melt  | 75.2                  | 248.1            |            |        |         |         |
| eis-melt  | 68.8                  | 158.2            |            |        |         |         |
| inhA-mut melt   |                       |                  |            |        |         |         |
| katG-mut melt   | 68.4                  | 184.6            |            |        |         |         |
| fabG1-mut melt  |                       |                  |            |        |         |         |
| ahpC-mut melt   |                       |                  |            |        |         |         |
| gyrA1-mutA melt   |                       |                  |            |        |         |         |
| gyrA1-mutB melt   | 72.3                  | 125.0            |            |        |         |         |
| gyrA1-mutC melt   |                       |                  |            |        |         |         |
| gyrA2-mutA melt   | 76.0                  | 207.9            |            |        |         |         |
| gyrA2-mutB melt   |                       |                  |            |        |         |         |
| gyrA3-mutA melt   |                       |                  |            |        |         |         |
| gyrA3-mutB melt   | 76.5                  | 128.0            |            |        |         |         |
| gyrA3-mutC melt   |                       |                  |            |        |         |         |
| gyrB2-mut melt  |                       |                  |            |        |         |         |
| rrs-mut melt  |                       |                  |            |        |         |         |
| eis-mutA melt   |                       |                  |            |        |         |         |
| eis-mutB melt   |                       |                  |            |        |         |         |

**Gambar 12. MTB TERDETEKSI (MTB DETECTED); Resistansi INH dan FLQ Rendah TERDETEKSI (DETECTED)**

| Test Result     | Analyte Result  | Detail           | Melt Peaks | Errors | History | Support |
|-----------------|---|------------------|------------|--------|---------|---------|
| Assay Name      | MTB-XDR   | Version          | 3          |        |         |         |
| Test Result     | <b>MTB DETECTED;</b><br><b>INH Resistance DETECTED;</b><br><b>FLQ Resistance DETECTED;</b><br><b>AMK Resistance DETECTED;</b><br><b>KAN Resistance DETECTED;</b><br><b>CAP Resistance NOT DETECTED;</b><br><b>ETH Resistance NOT DETECTED</b> |                  |            |        |         |         |
| Analyte Name    | Melt Peak Temperature   | Melt Peak Height |            |        |         |         |
| inhA-melt       | 76.6  | 278.9            |            |        |         |         |
| katG-melt       |   |                  |            |        |         |         |
| fabG1-melt      | 71.7  | 226.6            |            |        |         |         |
| ahpC-melt       | 69.0  | 42.9             |            |        |         |         |
| gyrA1-melt      |   |                  |            |        |         |         |
| gyrA2-melt      |   |                  |            |        |         |         |
| gyrA3-melt      |   |                  |            |        |         |         |
| gyrB2-melt      | 69.8  | 68.7             |            |        |         |         |
| rrs-melt        | 75.3  | 198.7            |            |        |         |         |
| eis-melt        |   |                  |            |        |         |         |
| inhA-mut melt   |   |                  |            |        |         |         |
| katG-mut melt   | 68.5  | 204.1            |            |        |         |         |
| fabG1-mut melt  |   |                  |            |        |         |         |
| ahpC-mut melt   |   |                  |            |        |         |         |
| gyrA1-mutA melt |   |                  |            |        |         |         |
| gyrA1-mutB melt | 72.9  | 88.0             |            |        |         |         |
| gyrA1-mutC melt |   |                  |            |        |         |         |
| gyrA2-mutA melt |   |                  |            |        |         |         |
| gyrA2-mutB melt |   |                  |            |        |         |         |
| gyrA3-mutA melt |   |                  |            |        |         |         |
| gyrA3-mutB melt |   |                  |            |        |         |         |
| gyrA3-mutC melt | 69.1  | 113.4            |            |        |         |         |
| gyrB2-mut melt  |   |                  |            |        |         |         |
| rrs-mut melt    |   |                  |            |        |         |         |
| eis-mutA melt   | 71.6  | 183.4            |            |        |         |         |
| eis-mutB melt   |   |                  |            |        |         |         |

Gambar 13. MTB TERDETEKSI (MTB DETECTED); Resistansi INH, FLQ, AMK, dan KAN TERDETEKSI (DETECTED)

The image shows two screenshots of the Xpert MTB/XDR software interface. The top screenshot displays the 'Test Result' tab with the following information:

- Assay Name: MTB-XDR
- Version: 3
- Test Result: **MTB NOT DETECTED**

The bottom screenshot shows a table with the following columns: Analyte Name, Melt Peak Temperature, and Melt Peak Height. The table lists various analytes, including wild-type and mutant strains.

| Analyte Name    | Melt Peak Temperature | Melt Peak Height |
|-----------------|-----------------------|------------------|
| inhA-melt       |                       |                  |
| katG-melt       |                       |                  |
| fabG1-melt      |                       |                  |
| ahpC-melt       |                       |                  |
| gyrA1-melt      |                       |                  |
| gyrA2-melt      |                       |                  |
| gyrA3-melt      |                       |                  |
| gyrB2-melt      |                       |                  |
| rrs-melt        |                       |                  |
| eis-melt        |                       |                  |
| inhA-mut melt   |                       |                  |
| katG-mut melt   |                       |                  |
| fabG1-mut melt  |                       |                  |
| ahpC-mut melt   |                       |                  |
| gyrA1-mutA melt |                       |                  |
| gyrA1-mutB melt |                       |                  |
| gyrA1-mutC melt |                       |                  |
| gyrA2-mutA melt |                       |                  |
| gyrA2-mutB melt |                       |                  |
| gyrA3-mutA melt |                       |                  |
| gyrA3-mutB melt |                       |                  |
| gyrA3-mutC melt |                       |                  |
| gyrB2-mut melt  |                       |                  |
| rrs-mut melt    |                       |                  |
| eis-mutA melt   |                       |                  |
| eis-mutB melt   |                       |                  |

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Tabel 17. Reaktivitas Analitis (Inklusivitas) untuk Galur Rentan Obat

| Sample: (Sampel)      | Keturunan Galur  | inhA         | katG         | fabG1        | oxyR- <i>ahpC</i> <sup>a</sup> | gyrA1        | gyrA2        | gyrA3        | gyrB2        | rrs          | eis          |
|-----------------------|------------------|--------------|--------------|--------------|--------------------------------|--------------|--------------|--------------|--------------|--------------|--------------|
| ( <i>M.bovis</i> BCG) | Tidak ditentukan | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | GAGAL (FAIL)                   | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) |
| <i>M.bovis</i>        | Tidak ditentukan | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | GAGAL (FAIL)                   | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) |
| MTB (AR2)             | 2                | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS)                   | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) |
| MTB (GD139)           | 3                | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS)                   | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) |
| MTB (AH1)             | 4                | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS)                   | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) |
| MTB (HR36)            | 5                | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS)                   | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) |
| MTB (HR37Rv)          | 4                | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | GAGAL (FAIL)                   | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) |
| <i>M.canetti</i>      | Tidak ditentukan | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | GAGAL (FAIL)                   | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) |
| <i>M.microti</i>      | Tidak ditentukan | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS)                   | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) | LULUS (PASS) |

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| Identitas Galur | Gen   | Mutasi yang Diharapkan    | MTB Terdeteksi (MTB Detected) | Tm Probe Mutan Terdeteksi (jumlah positif/diuji)           | Hasil RESISTANSI TERDETEKSI (RESISTANCE DETECTED) yang Benar (jumlah positif/diuji) |
|-----------------|-------|---------------------------|-------------------------------|--|---|
| 14-14194        | gyrA  | GAC 94 GCC                | 3 / 3                         | gyrA1-MutA,<br>gyrA2-MutA                                  | FLQ [3/3]   |
|                 | katG  | AGC 315 ACC               |                               | katG-Mut (3/3)   | INH [3/3]   |
|                 | inhA  | C -15 T                   |                               | inhA-Mut (3/3)   | INH, ETH [3/3]  |
| 15-14175        | katG  | AGC 315 ACC               | 3 / 3                         | katG-Mut (3/3)   | INH [3/3]   |
|                 | eis   | -10G/A                    |                               | eis-Mut (3/3)  | KAN [3/3]   |
| 15-14191        | katG  | AGC 315 ACC               | 3 / 3                         | katG-Mut (3/3)   | INH [3/3]   |
|                 | eis   | -10G/A                    |                               | eis-Mut (3/3)  | KAN [3/3]   |
| 16-05612        | katG  | AGC 315 ACC               | 3 / 3                         | katG-Mut (3/3)   | INH [3/3]   |
|                 | inhA  | C -15 T                   |                               | inhA-Mut (3/3)   | INH, ETH [3/3]  |
|                 | eis   | -12C/T                    |                               | eis-Mut (3/3)  | KAN [3/3]   |
| 16-05613        | katG  | AGC 315 ACC               | 3 / 3                         | katG-Mut (3/3)   | INH [3/3]   |
|                 | inhA  | C -15 T                   |                               | inhA-Mut (3/3)   | INH, ETH [3/3]  |
|                 | eis   | -12C/T                    |                               | eis-Mut (3/3)  | KAN [3/3]   |
| 14-13764        | katG  | AGC 315 ACC               | 3 / 3                         | katG-Mut (3/3)   | INH [3/3]   |
|                 | ahpC  | -48G/A                    |                               | ahpC-Mut (3/3)   | INH [3/3]   |
| 14-13806        | katG  | AGC 315 ACC               | 3 / 3                         | katG-Mut (3/3)   | INH [3/3]   |
|                 | ahpC  | -48G/A                    |                               | ahpC-Mut (3/3)   | INH [3/3]   |
| Klinis          | gyrA  | GCG 90 GTG,<br>GAC 94 GGC | 3 / 3                         | gyrA3-MutB (3/3)   | FLQ [3/3]   |
|                 | inhA  | C -15 T                   |                               | inhA-Mut (3/3)   | INH [3/3]   |
|                 | ahpC  | G-6A                      |                               | ahpC (2/3) <sup>b</sup>                                    | INH [3/3]   |
| Klinis          | katG  | AGC 315 ACC               | 3 / 3                         | katG Mut (3/3)   | INH [3/3]   |
| Klinis          | gyrB2 | ACC 539 AAC               | 3 / 3                         | gyrB2 WT <sup>c</sup>                                      | *Resistansi tidak terdeteksi [0/3]  |
|                 | rrs   | A1401G                    |                               | rrs-Mut (3/3)  | AMK, CAP, KAN [3/3]   |
|                 | gyrA  | GCG 90 GTG                |                               | gyrA1 MuB (3/3),<br>gyrA2 MutA (3/3),<br>gyrA3 MutB (3/3)  | FLQ [3/3]   |
|                 | ahpC  | g -6 a                    |                               | ahpC Mut (3/3)   | INH [3/3]   |
|                 | inhA  | C -15 T                   |                               | inhA Mut (3/3)   | INH, ETH [3/3]  |
| Klinis          | gyrA  | TCG 91 CCG                | 3 / 3                         | gyrA1-MutB (3/3),<br>gyrA2-MutA (3/3),<br>gyrA3-MutC (3/3) | FLQ [3/3]   |

| Identitas Galur | Gen  | Mutasi yang Diharapkan | MTB Terdeteksi (MTB Detected) | Tm Probe Mutan Terdeteksi (jumlah positif/diuji) | Hasil RESISTANSI TERDETEKSI (RESISTANCE DETECTED) yang Benar (jumlah positif/diuji) |
|-----------------|------|------------------------|-------------------------------|--|---|
|                 | inhA | C -15 T                |                               | inhA Mut (3/3)                                   | INH, ETH [3/3]  |

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## 21 Referensi

1. WHO. Global Tuberculosis Report 2019. [www.who.int/tb/publications/global\\_report](http://www.who.int/tb/publications/global_report)
2. WHO. 2018, Rapid Communication: Key changes to treatment of multidrug- and rifampicin-resistant tuberculosis.
3. WHO, 2019, Rapid Communication: Key changes to the treatment of drug-resistant tuberculosis
4. Sulis G, Pai M (2020) Isoniazid-resistant tuberculosis: A problem we can no longer ignore. *PLoS Med* 17(1): e1003023
5. Soumitesh Chakravorty, Sandy S. Roh, et al Department of Medicine, New Jersey Medical School, Newark, New Jersey, USA Detection of Isoniazid-, Fluoroquinolone-, Amikacin-, and Kanamycin-Resistant Tuberculosis in an Automated, Multiplexed 10-Color Assay Suitable for Point-of-Care Use *J Clin Microbiol* Jan 2017 Volume 55 Issue 1.
6. Centers for Disease Control and Prevention. Biosafety in microbiological and biomedical laboratories. Chosewood, LC dan Wilson, DE (eds) (2009). Nomor Publikasi HHS (CDC) 21-1112.
7. Clinical and Laboratory Standards Institute (sebelumnya National Committee for Clinical Laboratory Standards). Protection of laboratory workers from occupationally acquired infections; Approved Guideline. Dokumen M29 (lihat edisi terbaru).
8. Clinical and Laboratory Standards Institute (sebelumnya National Committee for Clinical Laboratory Standards). Laboratory Detection and Identification of Mycobacteria; Document M48A (lihat edisi terbaru).
9. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on the classification labeling and packaging of substances and mixtures amending and repealing, List of Precautionary Statements, Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

10. Occupational Safety and Health Standards, Hazard Communication, Toxic and Hazard Substances (March 26, 2012) (29 C.F.R., pt. 1910, subpt. Z).
11. Kent PT, Kubica GP 1985. Public Health Mycobacteriology—A Guide for Level III Laboratory, Centers of Disease Control, Atlanta, Publication no. PB 86-216546
12. Camus JC, et al. 2002. Re-annotation of the genome sequence of *Mycobacterium tuberculosis* H37Rv. *Microbiology* 148:2967–73.)
13. Zhang et al. 1996. Molecular basis for the exquisite sensitivity of *Mycobacterium tuberculosis* to isoniazid. *Proc Natl Acad Sci U S A.* 1996;93(23):13212–13216
14. Maruri et al. 2012. A systematic review of gyrase mutations associated with fluoroquinolone-resistant *Mycobacterium tuberculosis* and a proposed gyrase numbering system. *J Antimicrob Chemother.*
15. Vadwai V, Ajbani K, Jose M, Vineeth VP, Nikam C, Deshmukh M, Shetty A, Soman R, Rodrigues C. Can inhA mutation predict ethionamide resistance? *Int J Tuberc Lung Dis.* 2013 Jan;17(1):129-30. doi: 10.5588/ijtld.12.0511. Epub 2012 Nov 8. PubMed PMID: 23146620
16. Dunn JJ, Starke JR, Revell PA. 2016. Laboratory diagnosis of *Mycobacterium tuberculosis* infection and disease in children. *J Clin Microbiol* 54:1434–1441. doi:10.1128/JCM.03043
17. Diagnostic Standards and Classification of Tuberculosis in Adults and Children, *Am J Respir Crit Care Med* Vol 161. pp 1376–1395, 2000.
18. Seifert M, Catanzaro D, Catanzaro A, Rodwell TC (2015) Genetic Mutations Associated with Isoniazid Resistance in *Mycobacterium tuberculosis*: A Systematic Review. *PLoS ONE* 10(3): e0119628. doi:10.1371/journal.pone.0119628
19. Farhat MR, Jacobson KR, Franke MF, Kaur D, Sloutsky A, Mitnick CD, Murray M. 2016. Gyrase mutations are associated with variable levels of fluoroquinolone resistance in *Mycobacterium tuberculosis*. *J Clin Microbiol* 54:727–733. doi:10.1128/JCM.02775-15.
20. Gagneux S, Small PM. Global phylogeography of *Mycobacterium tuberculosis* and implications for tuberculosis product development. *Lancet Infect Dis.* 2007 May;7(5):328-37.
21. Banada, P. et.al. Containment of Bioaerosol Infection Risk by the Xpert MTB/RIF Assay and Its Applicability to Point of Care Settings. *Journal of Clinical Microbiology.* 2010.48:10. 3551-3557.

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© 2020–2023 Cepheid.
Cepheid AB
Cepheid Switzerland GmbH