Xpert[®] MTB/XDR

REF GXMTB/XDR-10

Upute za upotrebu

CE **IVD**

Izjave o zaštitnom znaku, patentima i vlasničkim pravima

Cepheid[®], the Cepheid logo, GeneXpert[®], and Xpert[®] are trademarks of Cepheid, registered in the U.S. and other countries. All other trademarks are the property of their respective owners.

THE PURCHASE OF THIS PRODUCT CONVEYS TO THE BUYER THE NON-TRANSFERABLE RIGHT TO USE IT IN ACCORDANCE WITH THESE INSTRUCTIONS FOR USE. NO OTHER RIGHTS ARE CONVEYED EXPRESSLY, BY IMPLICATION OR BY ESTOPPEL. FURTHERMORE, NO RIGHTS FOR RESALE ARE CONFERRED WITH THE PURCHASE OF THIS PRODUCT.

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Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name MTB-XDR Version 3						
Test Result	MTB DETECTED; INH Resistance NOT DETECTED; FLQ Resistance NOT DETECTED; AMK Resistance NOT DETECTED; KAN Resistance NOT DETECTED; CAP Resistance NOT DETECTED; ETH Resistance NOT DETECTED					
Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Analyte Name	Melt Peak Temperature	Melt Peak Height				
inhA-melt	76.3	292.5				
katG-melt	73.8	107.0				
fabG1-melt	71.5	242.0				
ahpC-melt	68.7	41.3				
gyrA1-melt	76.2	73.9				
gyrA2-melt	70.4	75.8				
gyrA3-melt	71.0	129.8				
gyrB2-melt	69.5	77.8				
rrs-melt	75.0	188.7				
eis-melt	68.5	145.3				
inhA-mut melt						
katG-mut melt						
fabG1-mut melt						
ahpC-mut melt						
gyrA1-mutA melt						
gyrA1-mutB melt						
gyrA1-mutC melt						
gyrA2-mutA melt						
gyrA2-mutB melt						
gyrA3-mutA melt						
gyrA3-mutB melt						
gyrA3-mutC melt						
gyrB2-mut melt						
rrs-mut melt						
eis-mutA melt						
eis-mutB melt						

Slika 8. MTB OTKRIVEN; otpornost na INH, FLQ, AMK, KAN, CAP i ETH NIJE OTKRIVENA (MTB DETECTED; INH, FLQ, AMK, KAN, CAP, and ETH Resistance NOT DETECTED)

Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name	MTB-XDR	Version	3			
Test Result	MTB DETECTED; INH Resistance DETECTED; FLQ Resistance DETECTED; AMK Resistance DETECTED; KAN Resistance DETECTED; CAP Resistance DETECTED; ETH Resistance DETECTED					
Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Analyte Name	Melt Peak Temperature	Melt Peak Height				
inhA-melt						
katG-melt						
fabG1-melt						
ahpC-melt						
gyrA1-melt	76.1	90.0				
gyrA2-melt	69.6	39.7				
gyrA3-melt						
gyrB2-melt						
rrs-melt						
eis-melt						
inhA-mut melt	70.9	259.6				
katG-mut melt	68.4	214.0				
fabG1-mut melt	75.9	181.1				
ahpC-mut melt	66.2	68.2				
gyrA1-mutA melt						
gyrA1-mutB melt						
gyrA1-mutC melt						
gyrA2-mutA melt						
gyrA2-mutB melt						
gyrA3-mutA melt						
gyrA3-mutB melt	76.0	125.0				
gyrA3-mutC melt						
gyrB2-mut melt	66.0	103.2				
rrs-mut melt	71.0	125.7				
eis-mutA melt	71.4	163.9				
eis-mutB melt						

Slika 9. MTB OTKRIVEN; otpornost na INH, FLQ, AMK, KAN, CAP i ETH OTKRIVENA
(MTB DETECTED; INH, FLQ, AMK, KAN, CAP, and ETH Resistance DETECTED)

Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name MTB-XDR Version 3						
Test Result	MTB DETECTED; INH Resistance DETECTED; FLQ Resistance NOT DETECTED; AMK Resistance NOT DETECTED; KAN Resistance NOT DETECTED; CAP Resistance NOT DETECTED; ETH Resistance NOT DETECTED					
Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Analyte Name	Melt Peak Temperature	Melt Peak Height				
inhA-melt	76.6	284.9				
katG-melt	74.0	105.2				
fabG1-melt						
ahpC-melt	69.0	35.4				
gyrA1-melt	76.6	65.2				
gyrA2-melt	70.4	64.9				
gyrA3-melt	71.4	92.2				
gyrB2-melt	69.7	84.7				
rrs-melt	75.3	146.8				
eis-melt	68.7	124.2				
inhA-mut melt						
katG-mut melt						
fabG1-mut melt	75.9	178.0				
ahpC-mut melt						
gyrA1-mutA melt						
gyrA1-mutB melt						
gyrA1-mutC melt						
gyrA2-mutA melt						
gyrA2-mutB melt						
gyrA3-mutA melt						
gyrA3-mutB melt						
gyrA3-mutC melt						
gyrB2-mut melt						
rrs-mut melt						
eis-mutA melt						
eis-mutB melt						

Slika 10. MTB OTKRIVEN; otpornost na INH OTKRIVENA (MTB DETECTED; INH Resistance DETECTED)

Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name	MTB-XDR	Version	4			
Test Result	MTB DETECTED; INH Resistance DETECTED; FLQ Resistance NOT DETECTED; AMK Resistance INDETERMINATE; KAN Resistance DETECTED; CAP Resistance INDETERMINATE; ETH Resistance DETECTED					
Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Analyte Name	Melt Peak Temperature	Melt Peak Height				
inhA-melt						
katG-melt						
fabG1-melt	71.5	254.6				
ahpC-melt	68.7	49.4				
gyrA1-melt	76.3	62.9				
gyrA2-melt	70.2	59.8				
gyrA3-melt	71.5	56.5				
gyrB2-melt	69.4	74.8				
rrs-melt						
eis-melt						
inhA-mut melt	70.9	277.7				
katG-mut melt	68.2	157.7				
fabG1-mut melt						
ahpC-mut melt						
gyrA1-mutA melt						
gyrA1-mutB melt						
gyrA1-mutC melt						
gyrA2-mutA melt						
gyrA2-mutB melt						
gyrA3-mutA melt						
gyrA3-mutB melt						
gyrA3-mutC melt						
gyrB2-mut melt						
rrs-mut melt						
eis-mutA melt						
eis-mutB melt	62.6	46.5				

Slika 11. MTB OTKRIVEN; otpornost na INH i KAN OTKRIVENA; NEODREĐEN REZULTAT za AMK i CAP (MTB DETECTED; INH and KAN Resistance DETECTED; AMK and CAP INDETERMINATE)

Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name MTB-XDR Version 3						
Test Result	MTB DETECTED; INH Resistance DETECTED; Low FLQ Resistance DETECTED; AMK Resistance NOT DETECTED; KAN Resistance NOT DETECTED; CAP Resistance NOT DETECTED; ETH Resistance NOT DETECTED					
Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Analyte Name	Melt Peak Temperature	Melt Peak Height				
inhA-melt	76.5	313.1				
katG-melt						
fabG1-melt	71.7	211.5				
ahpC-melt	69.0	47.2				
gyrA1-melt						
gyrA2-melt						
gyrA3-melt						
gyrB2-melt	69.6	81.1				
rrs-melt	75.2	248.1				
eis-melt	68.8	158.2				
inhA-mut melt						
katG-mut melt	68.4	184.6				
fabG1-mut melt						
ahpC-mut melt						
gyrA1-mutA melt						
gyrA1-mutB melt	72.3	125.0				
gyrA1-mutC melt						
gyrA2-mutA melt	76.0	207.9				
gyrA2-mutB melt						
gyrA3-mutA melt						
gyrA3-mutB melt	76.5	128.0				
gyrA3-mutC melt						
gyrB2-mut melt						
rrs-mut melt						
eis-mutA melt						
eis-mutB melt						

Slika 12. OTKRIVENA NISKA RAZINA MTB-a; OTKRIVENA otpornost na INH i niska otpornost na FLQ (MTB DETECTED; INH and Low FLQ Resistance DETECTED)

Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name MTB-XDR Version 3						
Test Result MTB DETECTED; INH Resistance DETECTED; FLQ Resistance DETECTED; AMK Resistance DETECTED; KAN Resistance DETECTED; CAP Resistance NOT DETECTED; ETH Resistance NOT DETECTED						
Analyte Name	Melt Peak Temperature	Melt Peak Height				
inhA-melt	76.6	278.9				
katG-melt						
fabG1-melt	71.7	226.6				
ahpC-melt	69.0	42.9				
gyrA1-melt						
gyrA2-melt						
gyrA3-melt						
gyrB2-melt	69.8	68.7				
rrs-melt	75.3	198.7				
eis-melt						
inhA-mut melt						
katG-mut melt	68.5	204.1				
fabG1-mut melt						
ahpC-mut melt						
gyrA1-mutA melt						
gyrA1-mutB melt	72.9	88.0				
gyrA1-mutC melt						
gyrA2-mutA melt						
gyrA2-mutB melt						
gyrA3-mutA melt						
gyrA3-mutB melt						
gyrA3-mutC melt	69.1	113.4				
gyrB2-mut melt						
rrs-mut melt						
eis-mutA melt	71.6	183.4				
eis-mutB melt						

Slika 13. OTKRIVEN MTB; OTKRIVENA otpornost na INH, FLQ, AMK i KAN (MTB DETECTED; INH, FLQ, AMK, and KAN Resistance DETECTED)

Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Assay Name	MTB-XDR	Version	3			
Test Result	MTB NOT DETECTED					

Test Result	Analyte Result	Detail	Melt Peaks	Errors	History	Support
Analyte Name	Melt Peak Temperature	Melt Peak Height				
inhA-melt						
katG-melt						
fabG1-melt						
ahpC-melt						
gyrA1-melt						
gyrA2-melt						
gyrA3-melt						
gyrB2-melt						
rrs-melt						
eis-melt						
inhA-mut melt						
katG-mut melt						
fabG1-mut melt						
ahpC-mut melt						
gyrA1-mutA melt						
gyrA1-mutB melt						
gyrA1-mutC melt						
gyrA2-mutA melt						
gyrA2-mutB melt						
gyrA3-mutA melt						
gyrA3-mutB melt						
gyrA3-mutC melt						
gyrB2-mut melt						
rrs-mut melt						
eis-mutA melt						
eis-mutB melt						

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Identifikacijska oznaka soja	Gen	Očekivana mutacija	MTB otkriven	Tm sonde mutanta otkriven (broj pozitivnih/testiranih)	Točni rezultati OTPORNOST OTKRIVENA (RESISTANCE DETECTED) (broj pozitivnih/testiranih)
Klinički	gyrA	GAC 94 TAC	3 / 3	gyrA1-MutB (3/3); gyrA3-MutC (3/3)	FLQ [3/3]
	katG	AGC 315 ACC		katG Mut (3/3)	INH [3/3]
	fabG1	G609A		fabG1 Mut (3/3)	INH [3/3]
Klinički	gyrA	GGC 88 GCC, GCG 90 GTG, TCG 91 CCG	3 / 3	gyrA1-MutB (2/3), ^a gyrA1-MutC (2/3), gyrA2-MutA (3/3), gyrA3-MutB (1/3)	FLQ [3/3]
	katG	AGC 315 ACC		katG Mut (3/3)	INH [3/3]
	rrs	A1401G		rrs-Mut (3/3)	AMK, CAP, KAN [3/3]
Klinički	gyrA	GAC 94 GGC	3 / 3	gyrA3-MutB (3/3)	FLQ [3/3]
	katG	AGC 315 ACC		katG Mut (3/3)	INH [3/3]
	rrs	A1401G		rrs-Mut (3/3)	AMK, CAP, KAN [3/3]
14-14194	gyrA	GAC 94 GCC	3 / 3	gyrA1-MutA, gyrA2-MutA	FLQ [3/3]
	katG	AGC 315 ACC		katG-Mut (3/3)	INH [3/3]
	inhA	C -15 T		inhA-Mut (3/3)	INH, ETH [3/3]
15-14175	katG	AGC 315 ACC	3 / 3	katG-Mut (3/3)	INH [3/3]
	eis	-10G/A		eis-Mut (3/3)	KAN [3/3]
15-14191	katG	AGC 315 ACC	3 / 3	katG-Mut (3/3)	INH [3/3]
	eis	-10G/A		eis-Mut (3/3)	KAN [3/3]
16-05612	katG	AGC 315 ACC	3 / 3	katG-Mut (3/3)	INH [3/3]
	inhA	C -15 T		inhA-Mut (3/3)	INH, ETH [3/3]
	eis	-12C/T		eis-Mut (3/3)	KAN [3/3]
16-05613	katG	AGC 315 ACC	3 / 3	katG-Mut (3/3)	INH [3/3]
	inhA	C -15 T		inhA-Mut (3/3)	INH, ETH [3/3]
	eis	-12C/T		eis-Mut (3/3)	KAN [3/3]
14-13764	katG	AGC 315 ACC	3 / 3	katG-Mut (3/3)	INH [3/3]

Identifikacijska oznaka soja	Gen	Očekivana mutacija	MTB otkriven	Tm sonde mutanta otkriven (broj pozitivnih/testiranih)	Točni rezultati OTPORNOST OTKRIVENA (RESISTANCE DETECTED) (broj pozitivnih/testiranih)
	ahpC	-48G/A		ahpc-Mut (3/3)	INH [3/3]
14-13806	katG	AGC 315 ACC	3 / 3	katG-Mut (3/3)	INH [3/3]
	ahpC	-48G/A		ahpc-Mut (3/3)	INH [3/3]
Klinički	gyrA	GCG 90 GTG, GAC 94 GGC	3 / 3	gyrA3-MutB (3/3)	FLQ [3/3]
	inhA	C -15 T		inhA-Mut (3/3)	INH [3/3]
	ahpC	G-6A		ahpC (2/3) ^b	INH [3/3]
Klinički	katG	AGC 315 ACC	3 / 3	katG Mut (3/3)	INH [3/3]
Klinički	gyrB2	ACC 539 AAC	3 / 3	gyrB2 WT ^c	*Nije otkrivena otpornost [0/3]
	rrs	A1401G		rrs-Mut (3/3)	AMK, CAP, KAN [3/3]
	gyrA	GCG 90 GTG		gyrA1 MuB (3/3), gyrA2 MutA (3/3), gyrA3 MutB (3/3)	FLQ [3/3]
	ahpC	g -6 a		ahpC Mut (3/3)	INH [3/3]
	inhA	C -15 T		inhA-Mut (3/3)	INH, ETH [3/3]
Klinički	gyrA	TCG 91 CCG	3 / 3	gyrA1-MutB (3/3), gyrA2-MutA (3/3), gyrA3-MutC (3/3)	FLQ [3/3]
	inhA	C -15 T		inhA-Mut (3/3)	INH, ETH [3/3]

- ^a Ovaj uzorak koji sadržava tri različite mutacije u genu *gyrA* nije dao mutantne Tm-ove za sonde sva tri gena *gyrA* svaki puta. No kako bi se točno odredila otpornost, najmanje jedna sonda mora dati mutantne Tm-ove. Rezultat je bio točan za sva ponavljanja budući da je najmanje jedna sonda za gen *gyrA* uvijek dala najmanje jedan Tm mutanta pri testiranju.
- ^b Ovaj je uzorak dvostruki mutant gena *katG* / *ahpC*. Ponavljanje s mutantnim Tm-om nedostajućeg gena *ahpC* nazvan je INH-R zbog prisutnosti mutacije *katG*, koja je otkrivena testom.
- ^c Ova specifična mutacija nije otkrivena testom. No postoje ograničeni klinički dokazi da ta mutacija zapravo može pridonijeti otpornosti fluorokinolona (FLQ) (mutacija s niskom pouzdanosti za otpornost na FLQ).

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21 Reference

1. WHO. Global Tuberculosis Report 2019 www.who.int/tb/publications/global_report
2. WHO. 2018, Rapid Communication: Key changes to treatment of multidrug- and rifampicin-resistant tuberculosis.
3. WHO, 2019, Rapid Communication: Key changes to the treatment of drug-resistant tuberculosis
4. Sulis G, Pai M (2020) Isoniazid-resistant tuberculosis: A problem we can no longer ignore. PLoS Med 17(1): e1003023
5. Soumitesh Chakravorty, Sandy S. Roh, et al Department of Medicine, New Jersey Medical School, Newark, New Jersey, USA Detection of Isoniazid-, Fluoroquinolone-, Amikacin-, and Kanamycin-Resistant Tuberculosis in an Automated, Multiplexed 10-Color Assay Suitable for Point-of-Care Use J Clin Microbiol Jan 2017 Volume 55 Issue 1.
6. Centers for Disease Control and Prevention. Biosafety in microbiological and biomedical laboratories. Chosewood, LC and Wilson, DE (eds) (2009). HHS Publication number (CDC) 21-1112.
7. Clinical and Laboratory Standards Institute (formerly National Committee for Clinical Laboratory Standards). Protection of laboratory workers from occupationally acquired infections; Approved Guideline. Dokument M29 (proučite najnovije izdanje).
8. Clinical and Laboratory Standards Institute (formerly National Committee for Clinical Laboratory Standards). Laboratory Detection and Identification of Mycobacteria; Document M48A (proučite najnovije izdanje).

9. UREDBA (EZ) br. 1272/2008 EUROPSKOG PARLAMENTA I VIJEĆA od 16. prosinca 2008. o razvrstavanju, označivanju i pakiranju tvari i smjesa, o izmjeni i stavljanju izvan snage Popisa oznaka obavijesti, Direktive 67/548/EEZ i Direktive 1999/45/EZ i o izmjeni Uredbe (EZ) br. 1907/2006.
10. Propisi o sigurnosti i zdravlju na radu, Priopćavanje opasnosti, Otrovne i opasne tvari (26. ožujka 2012.) (29 C.F.R., str. 1910, podtočka Z.
11. Kent PT, Kubica GP 1985. Public Health Mycobacteriology—A Guide for Level III Laboratory, Centers of Disease Control, Atlanta, Publication no. PB 86-216546
12. Camus JC, et al. 2002. Re-annotation of the genome sequence of *Mycobacterium tuberculosis* H37Rv. *Microbiology* 148:2967–73.)
13. Zhang et al. 1996. Molecular basis for the exquisite sensitivity of *Mycobacterium tuberculosis* to isoniazid. *Proc Natl Acad Sci U S A*. 1996;93(23):13212–13216
14. Maruri et al. 2012. A systematic review of gyrase mutations associated with fluoroquinolone-resistant *Mycobacterium tuberculosis* and a proposed gyrase numbering system. *J Antimicrob Chemother*.
15. Vadwai V, Ajbani K, Jose M, Vineeth VP, Nikam C, Deshmukh M, Shetty A, Soman R, Rodrigues C. Can inhA mutation predict ethionamide resistance? *Int J Tuberc Lung Dis*. 2013 Jan;17(1):129-30. doi: 10.5588/ijtld.12.0511. Epub 2012 Nov 8. PubMed PMID: 23146620
16. Dunn JJ, Starke JR, Revell PA. 2016. Laboratory diagnosis of *Mycobacterium tuberculosis* infection and disease in children. *J Clin Microbiol* 54:1434–1441. doi:10.1128/JCM.03043
17. Diagnostic Standards and Classification of Tuberculosis in Adults and Children, *Am J Respir Crit Care Med* Vol 161 pp 1376–1395, 2000.
18. Seifert M, Catanzaro D, Catanzaro A, Rodwell TC (2015) Genetic Mutations Associated with Isoniazid Resistance in *Mycobacterium tuberculosis*: A Systematic Review. *PLoS ONE* 10(3): e0119628. doi:10.1371/journal.pone.0119628
19. Farhat MR, Jacobson KR, Franke MF, Kaur D, Sloutsky A, Mitnick CD, Murray M. 2016. Gyrase mutations are associated with variable levels of fluoroquinolone resistance in *Mycobacterium tuberculosis*. *J Clin Microbiol* 54:727–733. doi:10.1128/JCM.02775-15.
20. Gagneux S, Small PM. Global phylogeography of *Mycobacterium tuberculosis* and implications for tuberculosis product development. *Lancet Infect Dis*. 2007 May;7(5):328-37.
21. Banada, P. et.al. Containment of Bioaerosol Infection Risk by the Xpert MTB/RIF Assay and Its Applicability to Point of Care Settings. *Journal of Clinical Microbiology*. 2010.48:10. 3551-3557.

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© 2020–2023 Cepheid.
Cepheid AB
Cepheid Switzerland GmbH