

# Cepheid C360

---



## Administrative Features Operator Manual



Nonmedical-device software  
system for disease surveillance

301-8332, Rev. J  
December 2023



# Legal Information

---

Cepheid<sup>®</sup>, the Cepheid logo, GeneXpert<sup>®</sup>, and Xpert<sup>®</sup> are trademarks of Cepheid, registered in the U.S. and other countries.

All other trademarks are the property of their respective owners.

© 2016-2023 Cepheid.

## Disclaimers

The Cepheid C360 software is a hosted platform that provides Internet-based services, including management, aggregation, sharing, alerting, and reporting data and real-time disease surveillance information (“**Cepheid C360 Software**”). Cepheid C360 Software enables customers to upload data from Cepheid compatible and connected systems, including the GeneXpert<sup>®</sup> family of systems (“**Cepheid Systems**”) in order to access this data, process it in various ways and share it with others.

The customer is responsible for ensuring that their data processing through the C360 Software complies with all applicable laws and regulations, including those concerning applicable privacy and data protection. The customer is also responsible for ensuring that such data processing is secure, including providing a secure site environment for and operation for the Cepheid System. The Cepheid C360 Software requires internet connectivity, and the customer is responsible for providing a secure internet connection and integrating the Cepheid C360 Software with the customer’s existing systems, networks and servers.

**To the fullest extent permissible by law, Cepheid excludes all warranties, representations, conditions or terms which may be implied. Cepheid does not warrant (i) the accuracy or completeness of any data accessed, tracked, collected or otherwise provided as part of the Cepheid C360 Software; or (ii) that the Cepheid C360 Software is free of bugs, speed issues or performance issues. The operation of the Cepheid C360 Software is dependent upon external factors including the proper and effective functioning of the internet, customer’s own and third party equipment and services, customer’s own configuration of the Cepheid C360 Software functionalities and customer’s own usage of the data made available. Cepheid gives no warranty regarding, and will have no responsibility or liability for, external factors.**

**C360 Sync is the software application on your instruments that established connection with Cepheid hosted servers.**

All examples (printouts, graphics, displays, screens, etc.) are for information and illustration purposes only and shall not be used for clinical or maintenance evaluations. Data shown in sample printouts and screens do not reflect actual patient names or test results. Labels depicted in the manual may appear different from actual product labels. Cepheid makes no representations or warranties about the accuracy and reliability of the information contained in this document. The information was developed to be used by persons trained and knowledgeable in the Cepheid Systems, including GeneXpert system operation or under the direct supervision of Cepheid technical support or service representatives. Updates to this document may be issued periodically.



Cepheid  
904 Caribbean Drive  
Sunnyvale, CA 94089  
USA

Telephone: + 1 408 541 4191  
Fax: + 1 408 541 4192

# Table of Contents

---

<b>Legal Information.....</b>	<b>3</b>
<b>1 Introduction.....</b>	<b>7</b>
1.1 About Cepheid C360 Administration.....	7
1.2 Related Publications.....	7
1.3 Technical Assistance.....	8
1.4 Addresses.....	8
<b>2 Features Overview.....</b>	<b>11</b>
2.1 Cepheid C360 Interface Overview.....	11
2.2 The Cepheid C360 Administrative Toolset Interface.....	12
2.3 Administrative Toolset Tables.....	13
2.4 Return to Previous Page.....	14
2.5 About Browser Variation.....	15
2.6 Cybersecurity.....	15
2.7 Logging In.....	15
2.8 Logging Out.....	18
2.9 Account Lockout.....	18
<b>3 Administering Cepheid C360.....</b>	<b>19</b>
3.1 Cepheid C360 Administration Setup.....	19
3.2 Verifying the Primary Institution Admin Account.....	19
3.3 Verify Institution Details.....	20
3.4 Adjusting Your Cepheid C360 Setup.....	22
3.4.1 Changing Basic User Settings.....	22
3.4.2 Changing Security Settings.....	23
3.5 Overview of the Cepheid C360 Configuration Settings.....	23
3.5.1 Setting Cepheid C360 Institution Policies.....	24
3.5.2 Setting Data-Collection Policies.....	25
3.5.3 Configure Instrument Registration Codes.....	27
3.6 About Updating ADFs and Instrument Software.....	29
3.7 Managing Labs.....	29
3.7.1 Adding a Lab.....	29
3.7.2 Verifying and Editing Lab Details.....	30
3.7.3 Deleting a Lab.....	31
3.7.4 Viewing a Lab's Assigned Instruments.....	31
3.7.5 Viewing a Lab's Assigned Users.....	32
3.8 Managing Users.....	32
3.8.1 Configuring Cepheid C360 User Roles.....	33
3.8.2 Adding Users.....	34
3.8.3 Adding Users In Batches.....	35
3.8.4 Editing User Information.....	40
3.8.5 Suspending or Deleting User Accounts.....	40

3.9 Overview of Cepheid C360 Audit Trails.....	41
3.9.1 Searching and Filtering Audit Trails.....	41
3.9.2 Viewing Audit Trail Details.....	43
<b>4 Supranational Institutions and Data-Sharing.....</b>	<b>45</b>
4.1 Making a Request To Share Data.....	45
4.2 Modifying a Sharing Request.....	49
4.3 Delete a Sharing Request.....	49
4.4 Responding to Sharing Requests.....	50
4.5 Stopping Data Sharing.....	51
4.6 Restarting Stopped Data Sharing.....	51

# 1 Introduction

---

This chapter introduces the *Cepheid C360 Administrative Features Operator Manual* and provides a brief description of the administrative features of Cepheid C360.

## 1.1 About Cepheid C360 Administration

Cepheid C360 is a web-based software application for administering Cepheid® systems and visualizing aggregated medical test data produced by Cepheid instruments.

This manual describes the following administrative functions of Cepheid C360.

- Administration of user accounts and passwords for logging in to Cepheid C360.
- Creation of laboratories.
- Assignment of user roles that restrict or permit institution personnel access to Cepheid C360 features depending on their job functions.
- Setting institutional policies that control a few overall settings for Cepheid C360
- Logging events taking place on all Cepheid instruments in the institution, including errors, routine maintenance tasks and other useful information.
- Maintaining information about instruments, laboratories and users that describes their relationships and allows patient test data to be recorded in a structured, traceable way.

## 1.2 Related Publications

For information about topics related to Cepheid C360, see the following publications.

<b>For this topic:</b>	<b>See this publication:</b>
The data-visualization features of Cepheid C360	<i>Cepheid C360 Data-Visualization Features Operator Manual</i>
Connecting instruments to Cepheid C360	<i>C360 Sync Installation and Networking Operator Manual</i>
Operating and administering the GeneXpert Xpress System	<i>GeneXpert Xpress System User's Guide</i>
Operating and administering GeneXpert Dx systems	<i>GeneXpert Dx System Operator Manual</i>
Operating and administering the GeneXpert Infinity System	<i>GeneXpert Infinity System Operator Manual</i>
Operating and administering the GeneXpert Edge System	<i>GeneXpert Edge System User's Guide</i>

### 1.3 Technical Assistance

Before contacting Cepheid Technical Support, collect the following information:

- Product name
- Instrument serial number or numbers
- Error messages (if any)
- Software version and, if applicable, Computer Service Tag number

#### United States Technical Support

Telephone: + 1 888 838 3222

Email: [techsupport@cepheid.com](mailto:techsupport@cepheid.com)

#### France Technical Support

Telephone: + 33 563 825 319

Email: [support@cepheideurope.com](mailto:support@cepheideurope.com)

Contact information for other Cepheid offices is available on our website at [www.cepheid.com/en/CustomerSupport](http://www.cepheid.com/en/CustomerSupport).

### 1.4 Addresses

#### Corporate Headquarters

Cepheid  
904 Caribbean Drive  
Sunnyvale, CA 94089  
USA

Telephone: + 1 408 541 4191

Fax: + 1 408 541 4192

[www.cepheid.com](http://www.cepheid.com)

#### European Headquarters

Cepheid Europe SAS  
Vira Solelh  
81470 Maurens-Scopont  
France

Telephone: + 33 563 825 300

Fax: + 33 563 825 301

[www.cepheidinternational.com](http://www.cepheidinternational.com)

#### Manufacturer Information





Cepheid  
904 Caribbean Drive  
Sunnyvale, CA 94089  
USA

Telephone: + 1 408 541 4191  
Fax: + 1 408 541 4192  
[www.cepheid.com](http://www.cepheid.com)



## 2 Features Overview

---

This chapter summarizes the purpose of Cepheid C360 and introduces the features and interface of the Cepheid C360 administrative toolset.

---

**Note** Cepheid C360 is a supplemental information-management web application that functions as an IT application. It is not involved in the diagnosis, treatment or monitoring of patients.

---

**Important** **Your institution is responsible for ensuring that the uploading of data to the Cepheid C360 server complies with applicable laws and regulations, including those concerning privacy and data-protection.**

---

In Cepheid C360, an institution is a single organization that can have a number of labs, institution administrators, laboratory administrators, and GeneXpert<sup>®</sup> Instrument Systems linked to their labs.

**Note** A supranational institution is an organization, such as the World Health Organization (WHO), that can send data-sharing requests to other institutions and, upon acceptance of such requests by the institutions, perform epidemiological interpretation of the aggregated data based on the data shared.


---

### 2.1 Cepheid C360 Interface Overview

In Cepheid C360, a column runs down the left edge of the screen. It displays the data-visualization tools near the top, and the administrative tools further down in the column.

The top-level items for data-visualization are **Medical, Instruments, Test Results** and **Notifications**. The data-visualization tools are explained in the *Cepheid C360 Data-Visualization Features Operator Manual*.


The Cepheid C360 administrative tools are located under the word **ADMIN** and are grouped into **Organization, Users, Logs, Audit Trails**, and **Configurations**. This operator manual explains the administrative tools and their menus.



## 2.2 The Cepheid C360 Administrative Toolset Interface

The **ADMIN** area includes five items, all of which open to provide search or configuration options. Some administrative items include several menus.

- The **Organization** item includes the **Institution** and **Laboratories** menus. They provide tools to configure and modify the institutions and laboratories in the organization. Note that the **Institution** menu displays only your institution. It is there for use by Cepheid Technical Support.
- **Users** is where administrators create or search for Cepheid C360 users, assign users to labs, disable or delete accounts, and perform other tasks affecting institution personnel.
- **Audit Trails** provides detailed information about what actions each Cepheid C360 user performs.
- The **Configurations** item includes **Registration Codes**, **Mobile Devices**, **Data Collection Policies**, and **Sharing Requests** menus. These tools are for an administrator to generate and manage registration codes that connect instrument systems to the Cepheid C360 system, configure what data Cepheid C360 collects from the instruments, and manage sharing requests between institutions.



Each of these tools is described in detail in this operator manual.

**Note**


Not all menu items are visible to all Cepheid C360 users. Access is based on [the role assigned to each user](#) in the Cepheid C360 system.

## 2.3 Administrative Toolset Tables

When using the Cepheid C360 administrative tools, the table view is generally where you access detailed information about anything and initiate necessary tasks. For example, under **Organization**, in the **Laboratories** table, you can search, filter and sort the list of labs, and export the list as a CSV file. Click an existing lab in the list to open its record and view or modify details about it, or click **New Laboratory** to create and configure a new lab. The search and filter tools available with each table are appropriate for the information the table displays. For example, the **Test Result Logs** table can be filtered by test GUID, instrument and result, among other settings, while the **Users** table can be filtered by lab, user role and related settings.

### Search and Filter Data

To search and filter data, look for the pull-down lists that generally appear above the column headers in the table itself. The search and filter menus vary, depending on the type of data in the table.



The screenshot shows a data table interface. At the top, there is a search bar and a download icon. Below the search bar, there are filter options for 'Institution' (set to 'CEP Medical Center 1 West'), 'Identifier' (with a 'Filter by identifier' link), and 'Name' (with a 'Filter by name' link). The table has four columns: 'Identifier', 'Name', 'Institution', and 'Address'. The 'Name' column header has a white arrow pointing up, indicating it is the current sort order. The 'Identifier' column header has a pale arrow pointing up, indicating it can be sorted. The table contains two rows of data. At the bottom, there is a pagination control showing 'Items per page: 25' and '1 - 2 of 2' with navigation arrows.

Identifier	Name ↑	Institution	Address
MED1-WEST-MOB	CEP Med 1 West Mobile Lab	CEP Medical Center 1 West	9680 W Burnside St Portland Oregon USA
MED1-WEST-STAT	CEP Med 1 West Stationary Lab	CEP Medical Center 1 West	9680 W Burnside St Portland Oregon USA

### Sort Data

Some tables' column headers are active, so that clicking one sorts the table by the data in that column. In the example shown below, the white arrow next to **Name** indicates that the table is currently sorted by name. The pale arrow next to **Email** appears when the mouse hovers over that column header and indicates that the table could be sorted on that value instead. Not all columns in all tables can sort data.

**Important** The search and filter menus for many Cepheid C360 tables include a field for selecting the institution, but institution personnel only have access to data from their own institution. The institution field is there for the use of Cepheid Technical Support.



The screenshot shows a table header with the following columns: Username, First Name, Last Name ↑, Email ↑, Phone Number, Institution, Laboratory, Role, and Status. The 'Last Name' and 'Email' headers have white arrows pointing up, indicating they are currently sorted. The 'Email' header also has a pale arrow pointing up, indicating it can be sorted.

Username	First Name	Last Name ↑	Email ↑	Phone Number	Institution	Laboratory	Role	Status
----------	------------	-------------	---------	--------------	-------------	------------	------	--------

### Download Data

To determine whether a table can be exported, look for the **Download CSV** icon, which is usually located near the top right corner of the screen.



Search, filter and sort the table to prepare the data you need. Then click the icon to generate a CSV file. On the **Test Results** screen in the Data-Visualization toolset, multiple CSV files can be stored in a cache until you are ready to download and save them. In all other cases the software opens a **Save** dialog immediately.

**Note** The CSV format is designed for import into spreadsheet, database and statistical software tools. In CSV files, commas separate each data point and paragraph breaks separate each row of data. Cepheid C360 can export up to 3,000 rows of a table as a CSV file. (The **Test Results** table in the data-visualization toolset can export up to 300,000 rows.)

Each section in this operator manual provides appropriate specifics about filtering, sorting and exporting tabular data.

## 2.4 Return to Previous Page

Navigating Cepheid C360 is straightforward. Some screens display the steps of your path through the interface, which you can retrace with a click.

---

---

[< Institutions](#) / [Institution Details](#) / Policies for CEP Medical Center 1 West

Because Cepheid C360 runs in a web browser, the browser's back arrow accomplishes the same thing and is available in all locations.

## 2.5 About Browser Variation

Because Cepheid C360 runs in a web browser, the exact appearance of the screen may differ slightly depending on which browser you use. However, the features and functions are identical in any browser.

---

**Important** Cepheid C360 is tested with the Chrome, Edge, and Firefox browsers. Cepheid does not test other browsers. Logging in to your Cepheid C360 system with another browser will display a warning message that the browser is not supported.

---

## 2.6 Cybersecurity

Cepheid C360 includes the following data-security and user-authentication safeguards:

- Data encryption during transport and at rest
- Strong password requirements
- Two-factor authentication
- Automatic logout after inactivity
- Timed lockout after failed logins
- Role-based access-control
- Audit logs that monitor user interactions

If you experience a cybersecurity event, contact Cepheid Technical Support.

## 2.7 Logging In

Follow these steps to log in to your Cepheid C360 account.

### First Login

When your account was created, you should have received a Welcome email with your user ID and a link to your institution's Cepheid C360 server. Follow the instructions in the email to complete setting up your account.

1. Follow the link in the Welcome email.

---

**Note** The Welcome email link is valid for only 24 hours.

---

2. On the **Welcome to C360** screen enter a password that meets the requirements displayed under **Password requirements**.

Welcome to C360

To complete your account setup, please enter a password and select two security questions

**Enter and confirm your new password**

Password requirements

- 14 –128 characters
- At least two lowercase letters
- At least two uppercase letters
- At least two numbers
- At least two special characters: (~!@#\$%^&\*\_-+.,/<>?:[]\|{})
- Cannot be your email or username address
- Cannot be one of your 10 previous passwords

New password \*

New password is required

Confirm password \*

**Select and answer your security questions**

Security question 1 \*

Security answer 1

Security question 2 \*

Security answer 2

Submit Cancel

3. Select two security questions under **Select and answer your security questions**. Provide answers for the questions that you will remember and that other people will not know.
4. Click **Submit**. You will be transferred to the regular login screen. Enter your user ID and new password, and your institution's alias. We recommend that you bookmark the Cepheid C360 site.

---

**Note**

The institution alias is provided by Cepheid to the institution's administration and is included in the Welcome email.

---

### Subsequent Logins

After your account is established, logging in to Cepheid C360 is similar to logging in to any system that uses two-factor authentication.

1. In a browser, open the URL for your institution's Cepheid C360 server.



2. Enter your username and the password you created when you completed setting up your account. Passwords are masked by default. To see the password as you enter it, click the eye icon next to the **Password** field. Passwords are case-sensitive, while usernames are not.
3. Enter the alias for your institution.

#### Note

If you enter an incorrect password or alias, your login fails but the error message doesn't explain why.

The login screen changes to display the **Code** field.

4. Look in your email or on your phone for a message from *Cepheid C360 team*. Enter the number provided in the **Code** field in the message and click **Login**.

### Replace Forgotten Password


If you forget your password you can replace it.

1. Click **Forgot password?** on the splash screen to open the **Forgot Password** screen.

2. Enter the email address associated with your Cepheid C360 account and click **Send Code**.
3. Look in your email for a message from *Cepheid C360 team*.
4. Enter the code provided in the message and click **Submit Code**.
5. [Create a new password](#). When changing your password during login, you do not set new security questions.

### 2.8 Logging Out

To log out of Cepheid C360, click the person icon next to your name, in the upper-right corner of the screen. When the menu displaying **Settings** and **Sign Out** opens, select **Sign Out**.



---

**Note** Cepheid C360 logs you out without a confirmation prompt.

---

Cepheid C360 logs you out automatically under the following circumstances:

- There has been no activity for 15 minutes.
  - The session has been active for more than 24 hours.
- Note**
- The browser window running the Cepheid C360 session is closed.
  - Your username and password are used to log in from another computer or from another browser on the same computer. This could happen if someone else logs in using your credentials, or if you log in from a second computer or browser.
- 

### 2.9 Account Lockout

If you make three unsuccessful attempts to log in within 30 minutes, the account locks and remains locked for 30 minutes.

## 3 Administering Cepheid C360

---

This chapter explains how to perform necessary administrative tasks in Cepheid C360. It covers the following tasks.

- How to create institutions and laboratories.
- How to create users and assign them to appropriate roles.
- How to set institution policies.
- How to view lab and instrument audit trails.
- How to list instruments, labs and personnel.
- How to change users' passwords, delete or suspend user accounts, and add users in batches.
- How to configure registration codes for connecting instruments to the C360 server and link them to their institution.

### 3.1 Cepheid C360 Administration Setup

During the sales process, a representative at your institution provided information to Cepheid, including the following:

- The name of your institution
- The names of one or more laboratories in the institution
- The name, email address and phone number of the user who will fill the primary Institution Admin role in Cepheid C360

To provision your institution's order, Cepheid performed the following tasks:

- Entered your institution's information into Cepheid C360
- Created login credentials for the primary institution administrator
- Sent a Welcome email with login information to the primary institution administrator

Upon receipt of the Welcome email, the primary institution administrator should log in to Cepheid C360, verify the information entered by Cepheid, add laboratories if needed, add users as needed, and set the institution policy.

---

**Note** The Welcome email link expires after 24 hours.

---

### 3.2 Verifying the Primary Institution Admin Account

When Cepheid hands off the institution it has created, the first thing you, as the primary institution administrator, should do is [ascertain that you can log in](#) and that your account in Cepheid C360 is set up correctly.

1. Under **ADMIN** click **Users** and find your name in the list. Under the **Role** column it should specify *Institution admin*.
2. Click your name.  
The **User Details** screen opens.

User	Institution
<b>Username</b> westadmin1	<b>Name</b> CEP Medical Center 1 West
<b>First name</b> CEPMedWest	<b>Location</b> Multnomah - Oregon - United States
<b>Last name</b> InstAdminOne	<b>Supranational Institution</b> No
<b>E-mail</b> instadmin1@cepmedwest.com	
<b>Phone</b>	
<b>Role</b> Institution admin	
<b>Created</b> 2/16/22, 5:15 PM	
<b>Last login</b> 3/6/22, 2:25 PM	

3. Review the information. If anything is incorrect or incomplete, click **Edit**.  
The **Edit Details** screen opens.
4. Correct the information, then click **Save**.

---

**Note**

You cannot change the user role or the institution's name or location or whether it is a supranational institution. If any of those fields are not correct, contact Cepheid Technical Support.

---

### 3.3 Verify Institution Details

After Cepheid, based on specifications provided by the institution administrator, creates an institution in Cepheid C360, the administrator should verify that it meets the specifications.

1. Click **Organization**, then click **Institution**.  
The **Institution Details** screen opens.

The screenshot displays the 'Institution' configuration page. At the top, the title 'Institution' is shown. Below it, the 'Mode' is set to 'Full' (indicated by a selected radio button). The 'Name' field contains 'Institution'. The 'Alias' field contains 'institution01'. The 'Cepheid Identifier' field also contains 'institution01'. The 'Display name' field contains 'Institution'. There is a checkbox for 'Supranational' which is currently unchecked. The 'Location' field contains 'Sweden'. The 'Assay code set' field is empty. The 'Resync Period (days)' is set to '2'. The 'Max Data Category Level' is set to 'Patient Identifier' (indicated by a selected radio button). At the bottom, there are three checked checkboxes: 'Show Module Info', 'Accepts sharing re...', and 'Enable Disease Su...'.

2. Review the name, alias, and supranational setting, as well as other controls on this screen such as whether to accept sharing requests. If Cepheid was responsible for setting up a lab at your institution, look for it in the **Laboratories** pane at the bottom of the **Institution Details** screen. You should also verify some other settings:
  - The resync period for your institution, that is, how frequently C360 Sync checks for test results that have not yet been uploaded to Cepheid C360.
  - Under **Max Data Category Level**, the level of detail about each test that the operator can configure data-collection policies to transfer from instruments to the Cepheid C360 database. The options range from simply whether the test returned a result, to what analytes it detected, to demographic data collected about the patient, or the patient ID.

### Important

**Your institution is responsible for ensuring that the uploading of data to the Cepheid C360 server complies with applicable laws and regulations, including those concerning privacy and data-protection. Discuss this with Cepheid before your institution is created in the Cepheid C360 system.**

- Whether to display information about instrument modules.
- Whether your institution can be solicited by supranational institutions to share data. **Accepts sharing requests** means that your institution is visible to supranational institutions and can receive share requests; it does not mean that your institution has agreed to share data.
- Disease surveillance and the **Medical Dashboard**. Disease surveillance must be active for Cepheid C360 to have access to the raw data returned by instruments before the GeneXpert software cleans up the test-results data.

**Note** If your institution was configured without enabling **Enable Disease Surveillance and Medical Dashboard**, the **Medical Dashboard** menu will not be visible. For details about disease surveillance and the **Medical Dashboard**, see the *Cepheid C360 Data-Visualization Features Operator Manual*.

---

**Note** If any of these settings are incorrect or you don't understand their impact, contact Cepheid Technical Support. Institution administrators cannot adjust settings on the Institution Details screen—only Cepheid Technical Support can make corrections.

---


## 3.4 Adjusting Your Cepheid C360 Setup

This section explains how to adjust a few basic settings of your Cepheid C360 account.


### 3.4.1 Changing Basic User Settings

Follow these steps to change details on the **User Settings** screen.

1. Click the person icon, then click **Settings**.  
The **User Settings** screen opens, with two buttons: **Change password** and **Edit**.



2. Click **Edit**.  
The email address and the phone number fields become editable.

A screenshot of the 'Institution Details' screen. At the top right, it shows 'InstAdmin User' and 'Institution admin' next to a person icon. Below this, there are two blue buttons: 'Change password' and 'Edit'. The main content area is titled 'Login information' and contains several fields: 'User' (Inst\_Admin\_role), 'E-mail' (instadmin4@cep.com), 'Phone', 'Alias' (MainInstitution01), 'Role' (Institution admin), and 'Last login' (5/15/23, 4:31 PM). Each field has a dotted line below it, indicating it is editable.

**Note** The only fields on the **User Settings** screen that can be edited after your account is set up are the email address and the phone number.

---

3. Make any required changes, then click **Save**. Click **Cancel** instead to revert to the previous email and phone.

### 3.4.2 Changing Security Settings

Follow these steps to change your account password and security questions.

1. Open the **User Settings** screen, then click **Change password**.  
The **Change Password** screen opens.

Please enter your current password to change your password and/or security questions.

Current Password

New Password

Set a new Password

Password requirements

- 14 –128 characters
- At least two lowercase letters
- At least two uppercase letters
- At least two numbers
- At least two special characters: (~!@#%&\*\_-=+,</>?;[]\|{})
- Cannot be your email or username address
- Cannot be one of your 10 previous passwords

New password

Confirm password

Security Questions

Set your security q...

Security question 1  
What was your childhood nickname? ▼

Security answer 1

Security question 2  
What is your oldest sibling's middle name? ▼

Security answer 2

2. Enter your current password.
3. Enter a new password, following the required password format shown on the screen.  
Confirm the new password on the second line.
4. Select or update one or both security questions from the pop-ups. Provide answers for the questions that you will remember and that other people will not know.
5. Click **Save**.


## 3.5 Overview of the Cepheid C360 Configuration Settings

This section describes how to configure general institution policies, set data-collection policies, and manage instrument-system registration codes.

### 3.5.1 Setting Cepheid C360 Institution Policies

The institution policy settings consist of a few general policies for how the institution is to operate. Only a user with the Institution Admin role can set institution policies.

1. Click **Organization**, then click **Institution**.  
The **Institution Details** screen opens.
2. Click **Institution Policy**.
3. In the **Institution Policy** screen, click **Edit**.
4. Use the table as a guide to make your policy choices.



Policy	Action
Automatically disable user account after days inactive?	Select to enter how many days an account can go without the user logging in before the account is disabled. The default value and maximum is 365 days.
Enable weaker password for lab assistant users?	This feature is no longer used.
Allow concealing a pathogen result from view?	Select to allow instrument operators to avoid displaying some results generated by an assay. This can prevent positive results of an endemic pathogen from overwhelming other positive results. (Not available for all assays or all instruments.)
Allow remote instrument software updates?	Select to allow software updates to be delivered over the internet. The default is selected.



Policy	Action
Bundle Download Code (BDC) required	<p>If <b>Allow remote instrument software updates?</b> is selected, the bundle download behavior is determined as follows:</p> <ul style="list-style-type: none"> <li>• If <b>checked</b> in Full mode, the Institution Admin receives the bundle code in the Download Center feature.</li> <li>• If <b>unchecked</b> in Full mode, no bundle code is sent.</li> <li>• If <b>checked</b> in SMO mode, the Business Administrator receives the bundle code in the Download Center feature.</li> <li>• If <b>unchecked</b> in SMO mode, no bundle code is sent.</li> </ul> <p>A bundle code can be used by an institution to download software or ADFs. The institution can determine if they require a bundle code to be entered prior to downloading software or ADFs. The bundle code expires in 48 hours.</p>

### 3.5.2 Setting Data-Collection Policies

Data-collection policies are useful to narrowly target specific information. For example, you can set up a data-collection policy to extract just the results of one assay conducted in a particular lab or on a particular instrument. Two data-collection policies could gather the same data from the same instruments but return different results based on different **Max Data Category Level** settings.

Data-collection policies are set up by Cepheid Technical Support or institution admins, usually at the request of an epidemiologist, to manage what information generated by instruments is presented in Cepheid C360 for data-visualization purposes. Supranational institutions do not receive data gathered by nonsupranational institutions' data-collection policies.

1. Under **Configurations**, click **Data Collection Policies**.

The **Data Collection Policies** screen opens, showing a list of your institution's currently active data-collection policies.

2. To create a new data-collection policy, click **New Data Collection Policy**.

The **New Data Collection Policy** screen opens. The **Data Collection Policy Enabled** slider defaults to being active. To create and implement a data-collection policy, leave the slider as is. You can set up and deploy a data-collection policy, then later make it temporarily inactive by shifting the slider to **Disabled**.

Data Collection Policy Details

Data Collection Policy Enabled


Name  
CEP West broad spectrum data

Start date  
3/6/2022

Institution  
CEP Medical Center 1 West

3. In the **Data Collection Policy Details** pane, give the policy a descriptive name and select the date when you want the policy to go into effect. The **Institution** field defaults to your institution.

When you complete entering information in the **Data Collection Policy Details** pane, the **Rules** pane opens just beneath it.



---

**Note** The **Patient Identifier** option may not be available for all customers as described in the provided user rights permissions document.

---

4. Fill out the fields in the **Rules** pane.
  - a) Select the level of data to collect.

The choices are constrained by the setting [in the Institution Details screen](#).
  - b) Check **Upload datagrams** to include datagrams with each test result.

---

**Note** Datagrams are a tool used by Cepheid Technical Support to monitor the integrity of data transmissions between instruments and the Cepheid C360 server. We recommend uploading datagrams except in situations of very limited bandwidth.

---

- c) Check **Sample ID** if you intend to collect that information from each test.

---

**Note** The **Sample ID** setting may not be visible, depending on how your institution was configured when first created by Cepheid.

---

- d) Click the line under **Assays** to open the list of assays. Choose **Select all** or click to select the assays you want to include. To add individual assays you must click the field line repeatedly to make the menu pop up. To remove an assay after you select it, click the small "X" next to the assay name. To remove them all, click **Select none**.
- e) From the **Fields** pop-up menu, choose **Select All**, or click to select one or more of the listed options. To select individual fields, click the field line repeatedly to make the menu pop up. To remove a field after you select it, click the small "X" next to the field name. To remove them all, click **Select none**.

Note: When you enter a patient's age, the **Age group** or **Age group CDC** (whichever is active based on the setup in C360 Sync) assigns the patient to standardized age categories for epidemiological purposes.
- f) Click **Add Rule** to append additional rules to the data-collection policy.

Rules

Max Data Category Level:  Result  Analytes  Patient Demographics  Patient Identifier

Upload datagrams  
If checked, datagram information will be included with each test, which provides additional information at the cost of increased bandwidth usage

Sample ID

Assays\*

Xpert C. difficile Xpert Flu/RSV XC (Flu only) Xpert Flu/RSV XC (Full) Select Assays

Select none

Fields

+ Add Rule

- Beneath the **Rules** pane, in the **Applies To** pane, narrow down the results by entering a location, one or more labs, or even specific instruments. The process is similar to adding assays to the data-collection policy.

Applies To

Select location

Multnomah - Oregon - United States Select Location Select none

Select laboratories

CEP Med 1 West Mobile Lab Select laboratories Select none

Select instruments

**Note** A data-collection policy can have as many rules as your institution has assays.

- Click **Save**.  
The new policy becomes operational and is added to the list of data-collection policies.
- To modify a data-collection policy, click the one of interest in the **Data Collection Policies** screen.  
The **Data Collection Policy Details** screen opens. You can change all of the policy's settings on this screen in the same way you created them.
- To delete a data-collection policy, open it from the **Data Collection Policies** screen. At the bottom of the **Data Collection Policy Details** screen click **Delete**.

### 3.5.3 Configure Instrument Registration Codes

Cepheid instruments require a registration code so that C360 Sync can connect the instruments to Cepheid C360. Follow these steps to generate and manage the registration codes. For more information on Cepheid C360 Sync see the *C360 Sync Installation and Networking Operator Manual*.

- Under the **Configurations** menu click **Registration Codes**.

The table of registration codes opens.

Name	Code	Institution	Applicable	Usage
	5cc04672-9e51-463d-8e09-338d2f2f7d31	CEP Medical Center 1 West	✓	4 / 50
West Med 1 Reg Code Aug-Dec	e5148bd6-18b6-4ff2-ba5d-2bcac6e11e85	CEP Medical Center 1 West	✓	0 / 50

Under the **Applicable** column in the table, registration codes displaying a white checkmark on a green background are active. Those displaying a white bar on a gray background are expired or deleted.

2. Click **New Registration Code** above the table of registration codes. The **New Registration Code** screen opens.

3. Give the new code a recognizable name. Set an expiration date in the **Expiration** date field and in the **Max usages** field specify how many Cepheid C360 Sync clients can use the registration code.

**Note** Hover over the information icons for details on the **Expiration** and **Max usages** fields.

**Note** The expiration and maximum usage settings are for security purposes.

4. Click **Generate**.  
The new code is created and appended to the end of the list in the **Registration Codes** table.

**Note** The **Registration Codes** table cannot be sorted, filtered or exported.

5. To modify a code's expiration date or usage limit, scroll through the list to find the code, then click to open it. Make the changes required, then click **Save**.

- To enter the registration code when registering a new instrument in Cepheid C360 Sync, you can copy it from the **Code** field to the Windows® Clipboard and then paste it into Cepheid C360 Sync.
- To delete an expired code, open it in the **Registration Codes** table, click **Delete**, then in the confirmation dialog click **Yes**.

## 3.6 About Updating ADFs and Instrument Software

Cepheid occasionally releases updated instrument software and new or updated assay definition files (ADFs). Based on your institution's discussions with Cepheid Technical Support, bundles of new software or updates will be prepared for your institution and delivered through C360 Sync, the software that manages communication between your instrument and the Cepheid cloud servers.

Downloading and installing ADFs and instrument software bundles from Cepheid C360 eliminates the need to receive and manage software distributed on CDs or DVDs. It allows Cepheid to ensure that your institution is using the optimal software for your instruments and workflow. Generally, it is the responsibility of an institution admin to download and install update bundles, though some of this task can be delegated to a laboratory admin. We recommend that the institution admin collaborate on this process with a system admin.

For detailed instructions on using C360 Sync to download and install bundles see the *C360 Sync Installation and Networking Operator Manual*.

## 3.7 Managing Labs

This section describes the tasks required to administer laboratories in Cepheid C360. The tasks include the following:

- Creating labs
- Assigning personnel and instruments
- Editing lab information
- Deleting labs

### Note

A lab in Cepheid C360 does not have to be a specific physical space or a division of the organization. You can create a lab in Cepheid C360 and assign instruments and personnel in whatever way you find convenient or necessary. A lab can be a temporary grouping, as needed to support a short-term project, or a permanent administrative unit. A lab can consist of instruments and personnel gathered in one room in a healthcare setting, or combined from different settings scattered across the globe.

### 3.7.1 Adding a Lab

Follow these steps to enter a lab at your institution into the Cepheid C360 system.

- Under **ADMIN**, click **Organization**, then click **Laboratories**.  
The **Laboratories** screen opens.
- Click **New Laboratory** near the top right corner of the screen.  
The **Add New Laboratory** screen opens.

< Laboratories / Add New Laboratory

Laboratory information

Name \*

Identifier \*

Type \*  
Stationary

Institution \*

Laboratory location

City \*

State \*

Country \*

Region \*

Location \*

Save Save and add another Cancel

3. Fill in the fields in the **Laboratory information** and **Laboratory location** panes.
  - For **Identifier**, define a brief label for the lab that colleagues will recognize. Lab identifiers can be up to 255 characters long, but we recommend using short, obvious identifiers.
  - Under **Type**, indicate whether the lab is mobile ("Ambulatory," such as a lab in a truck) or located in a healthcare facility ("Stationary").

---

**Note**

The Location field defaults to the same location as the institution that the lab belongs to. In the United States, this is always a county. Different countries have different geographical subdivisions.

---

4. Click **Save** or **Save and add another**.

### 3.7.2 Verifying and Editing Lab Details

Follow these steps to view and modify the settings for a laboratory.

1. Under **ADMIN**, click **Organization**, then click **Laboratories**.  
The **Laboratories** screen opens, displaying your institution's labs in a list. The list can be sorted by identifier, lab name, or institution.
2. Click the lab of interest.  
The **Laboratory Details** screen opens. It displays the same fields you filled out [when creating the lab](#).
3. Review the information on the **Laboratory Details** screen. Beneath the pane with information about the lab itself, there is another pane, called **Users**, with information about personnel assigned to the lab. The **Users** pane is closed by default. To open the list, click on it.

**Note** You cannot edit personnel in the Users list on this screen or sort the table. However, you can [export the list as a CSV file](#) by clicking the **Download CSV** icon just above the column headings.

4. If anything needs to be changed, click **Edit Laboratory**. Most fields on the **Edit Laboratory** screen can be modified. The institution and whether the lab is stationary or ambulatory cannot be modified, nor can the location once the lab has uploaded any data to the Cepheid C360 server.
5. Click **Update** to save the changes, or **Cancel** to revert to the previous settings. Click **Delete** to delete the lab, but [see additional considerations before deleting a lab](#).

### 3.7.3 Deleting a Lab

To delete a lab, select it from the **Laboratories** list, click **Edit Laboratory** in the **Laboratory Details** view, then click **Delete** at the bottom of the screen. When the confirmation dialog opens, click **Delete** again.

**A lab cannot be deleted under any of the following circumstances:**

- It has personnel associated with it.
  - It has instruments assigned to it.
- Important** • It has produced test results and uploaded the data to the Cepheid C360 database.

**You must remove personnel from the lab and ask Cepheid Technical Support to remove instruments before attempting to delete the lab. Only Cepheid Technical Support can remove instruments.**


Once a lab has produced test results and uploaded the data to Cepheid C360 the lab cannot be deleted, but it is possible to remove the lab's personnel and instruments. Uploaded data remains in Cepheid C360 after the instruments that produced it have been deleted.

### 3.7.4 Viewing a Lab's Assigned Instruments


Find the instruments assigned to a specific lab by following these steps.

**Note** The **Review Instruments** button on the **Instruments** screen opens a list of newly installed or reinstalled instruments waiting to be approved. That is a separate function from searching for a lab's operational instruments. For more information on installing instruments, see the *C360 Sync Installation and Networking Operator Manual*.

1. Click **Instruments** near the top of the **Data-Visualization** tools in the column at the left side of the screen.  
The **Instruments** screen opens.
2. From the **Model** pop-up menu at the top of the screen, select the instrument model or models, if desired.
3. Click the **Institution/Laboratory** pop-up. Check the checkbox next to your institution name.  
The right-pointing arrow points down and the institution's labs appear below the institution name. Each lab name is preceded by a checkbox.
4. Check the checkbox next to the lab of interest.  
Instruments assigned to that lab are listed in the table below the pop-up menus.



5. With the **Connection** pop-up you can further narrow down the list by searching for instruments based on how much time has elapsed since they last connected to Cepheid C360.



### 3.7.5 Viewing a Lab's Assigned Users

Follow these steps to generate a list of users assigned to a lab.

1. Under **ADMIN**, click **Organization**, then click **Laboratories**.  
The **Laboratories** screen opens with the table in the lower pane listing the laboratories at the institution.
2. Click the laboratory of interest.  
The **Laboratory Details** screen opens.
3. Scroll down and click the **Users** pane to display a table listing people assigned to the lab.  
The list of people cannot be sorted, and clicking a name does not open a detail screen about the user. However, the list can be searched, and it can be exported to a CSV file by clicking the **Download CSV** icon above the table.

## 3.8 Managing Users

This section describes tasks required to manage user accounts. This includes the following tasks:

- Adding, deleting and disabling user accounts
- Changing, deleting and resetting passwords
- Editing user contact information
- Assigning user roles



### 3.8.1 Configuring Cepheid C360 User Roles

Before adding users to Cepheid C360, you must first understand user roles.

Every Cepheid C360 user is assigned a user account with a user role. User roles define access privileges, just like the access privileges you may have on your organization's computer network.

Access privileges control what parts of the Cepheid C360 interface users can see as well as the data they can access and the tasks they can perform. Limiting access to the required features and data has the following advantages:

- It simplifies the interface for some users.
- Users cannot accidentally change someone else's work.
- Users cannot see information they are not authorized to see or perform tasks they are not authorized to perform.

Cepheid C360 has the following roles:

- **Institution Admin:** A user with the Institution Admin role has the same privileges as a laboratory admin. In addition, an institution admin can edit, add and delete labs, add additional institution admin and laboratory admin users, and set institution policy.
- **Laboratory Admin:** Users in the Laboratory Admin role can perform all administrative tasks for their own lab, but cannot add, delete or update other labs, add users with the Institution Admin role, or set institution policy.
- **Epidemiologist:** People in the Epidemiologist role view test results in the Cepheid C360 data-visualization interface and can change their own passwords. Epidemiologists do not run patient tests and cannot access or change any of the Cepheid C360 administrative features.

---

**Note** Supranational institutions have slightly different user roles, as described later in this manual.

---

#### User Role Privileges

The following table summarizes the tasks that different user roles can perform on the system.

Task	Institution Admin	Laboratory Admin	Epidemiologist
Log in to Cepheid C360	Yes	Yes	Yes
View Medical Dashboard	Yes	Yes	Yes
Add, invite edit, disable users	Yes	Yes (own lab)	No
Edit own password and security questions	Yes	Yes	Yes
Restore own forgotten password	Yes	Yes	Yes
Add, delete labs	Yes	No	No
Edit labs	Yes	Yes (own)	No
Edit institution policy	Yes	No	No
View institution policy	Yes	Yes	No
View lab installation report	Yes	Yes (own)	No
View Cepheid C360 audit	Yes	Yes (own)	No

Task	Institution Admin	Laboratory Admin	Epidemiologist
Download C360 Sync client	Yes	Yes	No
Create, edit, delete data-collection policy	Yes	No	No
Create, edit, delete registration codes	Yes	No	No
View sharing requests	Yes	Yes	No
Accept, reject, stop sharing requests	Yes	No	No
Create, edit, delete notifications	Yes	Yes	Yes
View notifications	Yes	Yes	Yes

### 3.8.2 Adding Users

Follow these steps to add a new employee to your Cepheid C360 environment.

1. Under **ADMIN** click **Users**.  
The **Users** screen opens.
2. Click **New User**, near the top right of the screen.  
The **New User** screen opens.
3. Under **User Details**, select a role for the user.  
Depending on the role selected, additional data-entry areas open in the **New User** screen.
4. Complete the additional fields as described below.
  - a) For a new person assigned the Laboratory Admin role, click the radio button for your institution under the **Filter by institution name** field under **User Access**. When you select the institution, the **Laboratory Assignment** box opens just to the right. Select the lab or labs for the new laboratory admin. Laboratory admins can be assigned to more than one lab in the list.

**Note** The displayed **Comparison levels** checkboxes vary with your selections.

- b) For a new person assigned the Institution Admin role, select the institution in the **User Access** area, as above. An institution can have more than one person with

the Institution Admin role. Institution admins have admin rights for all labs and are not assigned to a specific lab.

- c) For a new person assigned the Epidemiologist role, select the institution under **User Access**, as above. Then in the **Access by** pop-up menu that opens just to the right, select Everything, Location, or Laboratory. If you select Location, a new **Locations** entry field appears just below the pop-up menu. Enter a few letters to select from the list of locations associated with your institution. When you add an epidemiologist and select Location in the **Access by** pop-up, a **Comparison levels** setting appears beneath the institution name, with checkboxes for smaller divisions of the location that can include **Country, Region, District** and **Subdistrict**. The location divisions available vary by country and region. Note that epidemiologists do not need to be assigned to specific labs. However, if you select Laboratory in the **Access by** pop-up, the list of the institution's labs opens and the epidemiologist can be associated with one or more labs.

#### Note

Instead of the Epidemiologist role, supranational institutions have a Supranational Epidemiologist role. A person in that role can access data from more than one institution, depending on the institution's data-sharing arrangements. The **Access by** pop-up does not include Laboratory for the Supranational Epidemiologist role.

5. Complete the remaining fields, then click **Save** to accept the new record, **Save and add another** to continue adding employees, or **Cancel** to stop the process.


### 3.8.3 Adding Users In Batches

You can add up to 1500 users to Cepheid C360 in a single batch instead of one at a time. To do so, you must create a comma-separated list of the new users with all the items of information that will be uploaded to their **User Details** screens. There are three ways to do this: [in a template you download](#) from Cepheid C360, [in a text file](#), or [in Excel](#). Before starting, we recommend that you [read the notes at the end of this topic](#).

#### Overview

1. Under **ADMIN** click **Users**.

The **Users** screen opens.



Username	First Name	Last Name	Email	Phone Number	Institution	Laboratory	Role	Status
testadmin@cephid.se	F	testadmin			Cepheid		Lab assistant	👤
sohnadmin@cephid.se	F	sohnadmin			Sohna	SohnaLab1 - 3 more labs	Lab assistant	👤
thach_lab_assistant1005	Thach	ThachLabAssistant1005			Thach Institution 01	Thach Lab 01	Lab assistant	👤
thach_lab_assistant0007	Thach	ThachLabAssistant0007			Thach Institution 01	Thach Lab 01	Lab assistant	👤
thach_lab_assistant_002	thach	labassistant_002	hoang		Thach Institution 01	Thach Lab 01 - 1 more labs	Lab assistant	👤
thach_assistant_003	thach	assistant three	hoang		Thach Institution 01	Thach Lab 01	Lab assistant	👤
testadmin@cephid.se	F	testadmin			Cepheid		Lab assistant	👤
sohnadmin@cephid.se	F	sohnadmin			Sohna	SohnaLab1 - 3 more labs	Lab assistant	👤
username_565	subadmin	subadmin	sthriga_chen35@gmail.com		Sohna	SohnaLab1	Laboratory admin	👤
username_603	subadmin	subadmin	abovm.thuyng@gmail.com		Sohna	SohnaLab1	Laboratory admin	👤

2. Click **Add Users**.

Cepheid C360 displays the **New Users CSV** screen. This is the screen for uploading user batch files. The **Instructions for Bulk user upload** and the **Guidelines for editing the sample CSV file** drop-downs provide details on the procedure and on how to format the upload file.

### Use the Template

To enter the user information into a template, click **Download Sample CSV**. This will download a file called *userexample.csv* to your computer. You can rename this file if you desire.

**Note**


The *userexample.csv* file can be opened in a text editor such as Notepad or in Excel. For a small number of records it is easier to prepare the upload file in a text editor, as described later in this topic. Using a hybrid of Excel and Notepad for different stages of editing the *userexample.csv* file is more complex but allows you to enter numerous records in one step. This section describes the hybrid Excel/Notepad process.

1. Open the *userexample.csv* file in Excel.
2. Convert the field headings and sample data into columns by following these steps.
  - a. Click in the column A header to select the column.
  - b. Under the **Data** menu tab click **Text to Columns** to open the wizard.

	A	B	C	D	E	F	G	H
1	First Name	Last Name	Username	Phone Number	Email	Role	Laboratories	Password
2	labAdminFirstName	labAdminLastName	labAdminUsername	123456789	labAdminExample@email.com	laboratory_admin	exampleLaboratory1;ExampleLaboratory2;ExampleLaboratory3	ExamplePassword!!
3	labAssistantFirstName		labAssistantUsername	123456789		lab_assistant	exampleLaboratory1;ExampleLaboratory2;ExampleLaboratory3	
4	epidemiologistFirstName	epidemiologistLastName	epidemiologistUsername		epidemiologistExample@email.com	institution_epidemiologist		
5								

- c. In step 1 of the **Text to Columns** wizard check that the **Delimited** radio button is selected, then click **Next**.

- d. In step 2 of the **Text to Columns** wizard, under **Delimiters** uncheck the **Tab** box and check the **Comma** box.
- e. Click **Finish**. Do not click **Next**.



Data items in the rows are split into separate cells.

	A	B	C	D	E	F	G	H
1	First Name	Last Name	Username	Phone Number	Email	Role	Laboratories	Password
2	User	01	01user		01user@medwest.com	laboratory_admin	CEP Med 1 West Stationary Lab	
3	User	02	02user		02user@medwest.com	lab_assistant	CEP Med 1 West Stationary Lab	
4	User	03	03user		03user@medwest.com	institution_epidemiologist	CEP Med 1 West Stationary Lab	
5								
6								
7	labAdminFirstName	labAdminLastName	labAdminUsername	123456789	labAdminExample@email.com	laboratory_admin	exampleLaboratory1;ExampleLaboratory2;ExampleLaboratory3	ExamplePassword!!
8	labAssistantFirstName	labAssistantLastName	labAssistantUsername	123456789	epidemiologistExample@email.com	lab_assistant	exampleLaboratory1;ExampleLaboratory2;ExampleLaboratory3	
9	epidemiologistFirstName	epidemiologistLastName	epidemiologistUsername	123456789	epidemiologistExample@email.com	institution_epidemiologist		
10								

3. Import the rest of the user records you plan to upload to Cepheid C360. You can do this by copying and pasting or by using Excel's **Get Data** item on the **Data** menu tab. Depending on the source and format of the data you import, this step may require assistance from your institution's IT staff.
4. Save and close the CSV file.

#### Note


Excel automatically converts the comma-delimited data to tab-delimited format. Because Cepheid C360 can import only comma-delimited data, you must perform a few extra steps to change the delimiter back to commas.

5. Open the CSV file in Notepad.
6. In Notepad you will notice that the individual data items are widely separated. This is because there are tabs between them instead of commas. Select the blank space containing a tab between any two items on the page and copy it to the Clipboard.
7. From Notepad's **Edit** menu select **Replace...**
8. In the **Replace** dialog paste the tab from the Clipboard into the **Find what:** field.


#### Note

Notepad does not support using a special character, such as ^t, for a tab, and hitting the Tab key will simply move the cursor to the next item in the dialog.

9. Enter a comma in the **Replace with:** field.



10. Click **Replace All**.
11. Save and close the file.



The user records are now in comma-delimited format and the file is ready to upload to Cepheid C360, as described later in this topic.

#### Create a Text File

To enter information in a text file, create a separate row (a paragraph) for each person. Separate each item of information in the row with a comma. Include the following items in the order shown.

- First name
- Last name
- Username
- Phone number
- Email address
- Role
- Lab or labs
- Password

#### Use Excel

1. Create a new Excel file and save it in CSV format.
2. Use a separate row for each person.
3. Enter each data point in a separate cell.

---

**Important** The Excel upload file must be encoded as UTF-16LE. Excel supports UTF-16LE on the Macintosh but not on the Windows platform. Therefore, we do not recommend attempting to use Excel for the bulk upload process on a Windows system.

---

## Upload the File

1. When the file is ready, select your institution from the **Institution** field.
2. Click **Choose file**.
3. Find the file in the dialog and click **Open**.
4. In the **Add Users** screen click **Save**.

The names and their associated details are added to the list of users.

Several error conditions can generate an error message when you attempt to upload the file. Examples include a lab name that doesn't match any lab at the institution and a username with too few characters.

The screenshot shows a web interface titled "New Users CSV". It includes a "Select users CSV file:" section with a "Choose File" button and "No file chosen" text. To the right is a "Download Sample CSV" button. Below this is an "Institution" dropdown menu. At the bottom, a red error message box is displayed with the text: "CSV file errors" and "username is either too short, too long or contains illegal characters Row: 3 - USER\_VALIDATION\_TOO\_SHORT\_TOO\_LONG\_OR\_CONTAINS\_ILLEGAL\_CHARACTERS".

## Notes

### Note

The first row in the *userexample.csv* file is the field names (First Name, Last Name and so on). Do not erase this row. If you create a text file include a paragraph at the top with the correct number of commas, or enter the field names. If you create an Excel file include a row of blank cells, or enter the field names. The bulk user import ignores content in the first row of the upload file, so if the row of field names is missing from the *userexample.csv* file, or the text or Excel file do not provide an extra row, the import process will skip the first name in the file.

### Note

The *userexample.csv* file also includes several rows of example data. You can delete these rows to avoid including them in the upload.

### Note

Roles must be entered exactly as shown below. We recommend that you copy and paste role names from the template file.

- laboratory\_admin
- institution\_epidemiologist

### Note

Laboratories must be spelled exactly as they appear in your institution. We recommend that you copy and paste lab names from the list of labs in Cepheid C360.

### Note

To add a new user to more than one lab, include multiple labs between the same set of commas, or in the same cell if you are using Excel, and separate the labs with semi-colons.

### Note

If an item of information is not available, include an extra comma where that item belongs, or leave a blank cell if you are using Excel. All items are mandatory except a phone number.

### Note

Enter phone numbers with a + sign and the country code.

#### 3.8.4 Editing User Information

Follow these steps to modify information about an existing Cepheid C360 user.

1. Open the **Users** screen.
2. In the **Search** area, enter a username, email address, or first or last name in the **Search** field. Select a role from the **Role** pop-up menu, and select whether you're searching for enabled or disabled employee accounts from the **Status** pop-up menu.

**Note**

Hover over the information icon for a reminder of what criteria can be entered in the **Search** field.

---

**Note**

You can search and filter on more than one criterion.

---

The employees who meet the search criteria are listed in a table below the search area.

3. Click the employee you're seeking.  
The **User Details** screen opens.
  4. Click **Edit**.
  5. In edit mode you can change all information on the **User Details** screen except the institution and location under the **User Access** area.
- 

**Note**

If a new user has received a Welcome email but has not logged in, admins will see a **Resend Invitation** button next to the **Edit**, **Disable** and **Delete** buttons. Click **Resend Invitation** to remind the new user to log in.

---

#### 3.8.5 Suspending or Deleting User Accounts

To temporarily or permanently remove a user from Cepheid C360 follow these steps.

1. Open the **Users** screen and search for the user of interest.
2. In the table of users retrieved by the search, click the name of the user you're looking for to open that person's **User Details** screen.
3. To temporarily suspend the user's account, click **Disable**.

**Note**

Disable temporarily suspends the user's account, but the account remains in the system and can be reactivated when desired. To reactivate a disabled user, simply return to the **User Details** screen and click **Enable**

---



4. To permanently remove the user from the system, click **Delete**.

**Note**

Delete permanently removes the user from Cepheid C360. To reinstate a deleted user, you must recreate the account, following the same steps as for a new user.

---



## 3.9 Overview of Cepheid C360 Audit Trails


This section describes Cepheid C360 audit trails. It explains their purpose and how to use them.

Audit trails record detailed information about the actions each Cepheid C360 user performs. The audit trail interface provides a rich set of search and filter tools.

### 3.9.1 Searching and Filtering Audit Trails

Follow these steps to retrieve audit information from the database of recorded events.

1. Click **Audit Trails** in the left side menu under **ADMIN**.
2. In the **Audit Trails** screen, before beginning a search, click **Clear filters**. Any previous search criteria are removed.



3. In the **Search** pane, enter criteria for the records you want to retrieve.

**Note** You can use any combination of the filters to narrow down the results.

- a) Make selections as needed from the **Laboratory** and **User** pop-up menus.
- b) Select an audit type from the **Entity Type** pop-up menu. There are about 30 entity types to choose from, including assay, data-collection policy, registration code, sharing request, and test result.
- c) Select a value for the entity type from the **Activity Type** pop-up menu. Values available in the **Activity Type** menu depend on the entity type selection. For example, when **Entity Type** is set to Assay, the only item under **Activity Type** is Download CSV. When **Entity Type** is set to Data Collection Policy, the **Activity Type** choices are Create, Update, and Delete. And when **Entity Type** is set to Sharing Request, the **Activity Type** choices are Create, Update, Delete, Accept, and Reject. See the table below for details.

**Note** There may be audit types that you don't recognize. Not all audit types apply to all instrument models.

Entity Type	Activity Types
Additional fields	Create, update, delete
Notification definition	Create, update, delete
Assay	Download CSV
Assay code sets	Create, update, delete
Assay group	Update

Entity Type	Activity Types
Audit	Download CSV
Bundle	Submit, publish, reject, disable, verify, download CSV
Bundle deployment rule	Create, update, delete, download CSV
Data collection policy	Create, update, delete
Document	Create, send document, delete
Institution	Create, update, delete, download CSV, update institution policy, create institution policy
Instrument	Create, update, delete, download CSV, refresh configuration item, view system configuration, view system files, view Sync event logs, view rejected updates, download file, force resend of test results, refresh policies, retrieve infinity logs, retrieve system logs, request file, view policies, force resync of policies, force resync of test results, accept, reject, approve
Laboratory	Create, update, delete, download CSV
Registration code	Create, update, delete
Sharing request	Create, update, delete, accept, reject
Tech collection project	Create, update, delete
Test result	Details viewed, visibility changed, Cepheid internal assay details viewed, Cepheid internal assay visibility changed, download test result archive file, download CSV, Cepheid internal assay download CSV, Cepheid internal assay download GXX
Test result log	Download CSV
Test results recurrent exports	Create, update, delete
User	Login, logout, logout inactive, logout session duration exceeded, login failure, adopted profile, authorization failure, password change, create, update, delete, enable, disable, forgot password, password reset, password reset by admin, download CSV, security question change, resend invitation
Self test data	Download CSV
Instrument log	Download CSV

- d) Enter an additional criterion, if desired, in the **Search entity name or entity id or client ip** field. It is a free-text search field that can filter by instrument models, user names, telephone numbers, locations, or any other information captured in the audit database.
- e) Set start and end dates.

4. Click **Search**.

The records that meet the filter criteria are listed in a table below the **Search** pane. Sort the list by clicking any of the table's column headers, including the type of event, when it occurred, or the user involved.

### 3.9.2 Viewing Audit Trail Details

You can review details in an audit trail and export a batch of audit trails by following these steps.

1. Click an individual audit in the table to open it in the **Audit Details** screen. This screen has two panes, the **Audit trail details** pane and the **Audit trail changed attributes** pane.

Information on the **Audit trail details** pane varies slightly with the type of audit. Most include the event type, the person involved, a time stamp, and similar basic information.

Audit trail details	
Timestamp :	Mar 7, 2022, 10:06:35 AM
Entity Type :	User
Activity Type :	Disable
Entity Name :	CEPMedWest EpidemiologistTwo
Entity Id :	westepid2
User :	westadmin1
Access Role :	Institution admin
User's Institution GUID :	7cfb7a5a-4969-4bac-b40e-76b3770ee94d
Institution :	CEP Medical Center 1 West
Institution GUID :	7cfb7a5a-4969-4bac-b40e-76b3770ee94d
Client IP :	100.21.142.78
Source :	C360

2. Scroll down to the bottom of the **Audit Details** screen to view the **Audit trail changed attributes** pane. This provides specific details about the audit you are reviewing. For example, for a data-collection policy audit item, the **Audit trail changed attributes** shows the old rule and the new rule, while for a login failure it shows what occurred to prevent the person from logging in.

Audit trail changed attributes		
Attribute	Previous Value	New Value
Name	CEP Medical Center 1	CEP Medical Center 1 West

### 3 Administering Cepheid C360

---

---

---

**Note** Not all audit types involve changed attributes. For those audit types the **Audit trail changed attributes** pane does not appear.

---

3. Click the **Download CSV** icon to export the filtered and sorted list to a CSV file.

---

**Note** The **Audit Trails** CSV file includes only the information in the table. Information from the **Audit Details** screens is not included.

---

## 4 Supranational Institutions and Data-Sharing

---

This chapter describes how supranational institutions set up and send data-sharing requests, and how nonsupranational institutions manage the requests.

An institution defined as supranational in Cepheid C360 is an institution that aggregates medical test data from other medical centers. It may be part of a public-health agency such as the CDC or the WHO, or an entity of a private medical network tasked with gathering and interpreting information from the network's other institutions. Supranational institutions do not operate their own labs. Instead, their administrators request data from other institutions, and their epidemiologists interpret the data.

---

**Note** An institution is defined as supranational when Cepheid creates it in Cepheid C360. Once an institution has been defined as supranational it cannot be converted to a regular institution.

---

**Note** An institution can be enrolled in Cepheid C360 twice, as a supranational institution and a nonsupranational institution. That would allow it to issue sharing requests and to accept sharing requests from other institutions. To configure two versions of an institution, provide different identifiers and slightly different names. Keep in mind that supranational institutions cannot have laboratories.

---

To gather medical test data, a supranational institution sends out sharing requests that specify what information it is seeking. A sharing request might specify, for example, the results of a set of assays for a particular date range. The receiving institution or institutions can accept or reject the request. However, only supranational institutions and Cepheid Technical Support can initiate data-sharing requests.

---

**Note** The Supranational Admin and Supranational Epidemiologist roles have slightly different privileges in Cepheid C360 than the Institution Admin and Epidemiologist roles have. Generally, the supranational roles do not provide access to features that are internal to an institution, such as adding and inviting new institution personnel, adding and editing labs, viewing instrument configurations, and setting institution policy. In the same way, supranational admins and epidemiologists have a few privileges that institution personnel do not, such as adding and inviting new supranational personnel, and issuing or modifying sharing requests.

---

Data collected with data-sharing requests becomes available to the Cepheid C360 data-visualization tools. However, the shared data provided to the supranational institution includes only the test results collected by the sharing institution. Lab information, personnel, and other information is not included.

---

**Note** A supranational institution cannot distribute data gathered with sharing requests to other institutions.

---

### 4.1 Making a Request To Share Data

Administrators and epidemiologists at supranational institutions can follow these steps to generate sharing requests. Administrators and epidemiologists at nonsupranational institutions do not generate sharing requests, but they should understand how sharing requests work.

## 4 Supranational Institutions and Data-Sharing

### Note

For a nonsupranational institution to participate in sharing data, it must be configured to accept sharing requests at the time the institution is created in Cepheid C360. An institution configured to accept data-sharing requests is visible to supranational institutions that may send sharing requests. The institution is not obligated to accept them; it can accept or reject individual sharing requests. An institution not configured to accept data-sharing requests does not appear in the list of potential sharing institutions that supranational institutions draw from.

1. Under **ADMIN**, click the arrow next to the **Configurations** item to expand it. Then click **Sharing Requests**.

The list of previously created sharing requests opens in a table. The table can be sorted by institution, or alphabetically by each request's name.

Name	Institution	Assays	Periodicity	Data category	Granularity	Data from
BCR 2020s overview	CEP Supranational Medical East	Xpert BCR-ABL Monitor, Xpert BCR-ABL Ultra, Xpert BCR-ABL Ultra p190	Monthly	None	Laboratory	1/1/20
C.Diff Metro Portland	CEP Supranational Medical East	Xpert C. difficile	Daily	None	Test Result with Patient Demographics	1/1/22
Carba-R West	CEP Supranational Medical East	Xpert Carba-R	Real time	None	Laboratory	1/1/21
Fil and FV urban areas	CEP Supranational Medical East	Xpert Fil & FV	Daily	None	Country	1/31/21
GBS south	CEP Supranational Medical East	Xpert GBS LB, Xpert GBS	Real time	None	Test Result with Patient Demographics	2/14/22
Hemo fever	CEP Supranational Medical East	Xpert Ebola	Daily	None	Laboratory	1/31/22
Multnomah schools	CEP Supranational Medical East	Xpert Norovirus	Daily	None	Test Result with Patient Demographics	1/1/21
Recent STD incidence	CEP Supranational Medical East	Xpert HPV 16_18-45, Xpert HPV HR, Xpert HIV Qual, Xpert HIV Quant	Daily	None	Subdistrict	1/1/22
Recent Throat and Lung	CEP Supranational Medical East	Xpert Flu/RSV XC (Flu only), Xpert Flu, Xpert MRSA/SA Nasal Complete, Xpert MRSA, Xpert Xpress SARS-CoV-2/Flu/RSV, Xpert Xpress Strep A	Daily	None	Test Result with Patient Demographics	1/31/22

Items per page: 25 1 - 9 of 9 < >

2. Click **Create new sharing request**.

The **New sharing request** screen opens.

### Note

Only Supranational Institution Admin users can create a sharing request.

< Sharing Requests / New sharing request CEPSupra InstAdminOne  
Supernational admin

**Properties**

Name \*  
Recent Throat and Lung

Request data since \*  
1/31/2022

Assays \*

Xpert Flu/RSV XC (Flu only) Xpert Flu Xpert MRSA/SA Nasal Complete Xpert MRSA  
Xpert Xpress SARS-CoV-2/Flu/RSV Xpert Xpress Strep A Select none

Assays that currently belong to other sharing request will not be available for selection.

**Applies to \***

Locations  Institutions Add institutions

Search: CEP Medical Center 1 West, CEP Medical Center 2 South

**Data aggregation options**

Granularity:  Country  Region  District  Subdistrict  Laboratory  Test result  Test result with patient demographics

Time aggregation:  Real time  Daily  Weekly  Monthly

Get data 1 | days after the end of each day

3. Fill out the fields in the **Properties** pane.
  - a) Enter a recognizable name for the sharing request in the **Name** field, such as Influenza incidence in metropolitan New York.
  - b) Click the calendar icon next to the **Request data since** field and select a start date.  
If the start date is in the past, the sharing request retrieves the data retroactively.
  - c) Either click in the **Assays** field and select the tests of interest from the menu of choices, or click **Select all** in the menu.  
If you select an assay by mistake, click the small "X" to its right to remove it from the request.

**Important** You can send more than one sharing request to the same institution and you can include more than one assay in a request, but each assay can only be included in one request. To retrieve test results for a particular assay from multiple sharing institutions, include them all on the request. We recommend that you plan sharing requests carefully to avoid duplication.

4. In the **Applies to** pane select either the **Locations** or the **Institutions** radio button.
  - a) Click in the Search field. A menu of locations or institutions opens, depending on which radio button you selected.
  - b) Click the checkboxes of the locations or institutions you want to add to the sharing request. You can select multiple locations or institutions. There is no limit to the number of institutions included in a request. If you select by location, only the institutions in that location and that are configured to accept sharing requests receive the request.

**Note** If some of the institutions in a sharing request decline the request, data will still be retrieved from those institutions that accept the request.

---

5. Set the following details in the **Data aggregation options** pane.
  - a) Choose one radio button for **Granularity**. The choices are as follows:
    - Country
    - Region
    - District
    - Subdistrict
    - Laboratory
    - Test result
    - Test result with patient demographics

**Note** If an institution does not collect patient demographics or is set up to prevent sharing that information, it will not be included in the data shared after the request is accepted even if the requesting institution includes that in the sharing request.

---

**Note** Not all geographical areas provide the same levels of granularity. For details about your location contact Cepheid Technical Support.

---

- b) Choose one radio button for **Time aggregation**. The choices are as follows:
    - Real time
    - Daily
    - Weekly
    - Monthly
  - c) The choices for the final setting, **Get data "X" days after...**, vary depending on the **Time aggregation** setting. For example, if you select the **Weekly** radio button and enter 2 in the **Get data "X" days after...** field, data from institutions that accept the sharing request will be retrieved two days after the end of each week. These settings are visible when the institution admin at the sharing institution reviews the request on the **Sharing Request Details** screen.
6. Click **Save**.

The sharing request is added to the table of requests at the supranational institution and at the institution or institutions receiving the request. However, it will not be active until the specified institution or institutions accept the request.

**Note** After a sharing request has been created, a new pane, **Requested acceptance status**, appears at the bottom of the screen. Click the arrow next to the title to see whether the listed institutions have accepted the request.

---



Requested acceptance status	
<b>SathiyalInstitution</b>	
Xpert Xpress SARS-CoV-2/Flu/RSV	PENDING
Xpert Xpress SARS-CoV-2/Flu/RSV	PENDING
<b>STG Institution</b>	
Xpert Xpress SARS-CoV-2/Flu/RSV	PENDING
Xpert Xpress SARS-CoV-2/Flu/RSV	PENDING
<b>Thach Institution 01</b>	
Xpert Xpress SARS-CoV-2/Flu/RSV	PENDING
Xpert Xpress SARS-CoV-2/Flu/RSV	PENDING

## 4.2 Modifying a Sharing Request

After a request has been created, some of the parameters can be changed and some cannot. To change a request, open it by simply clicking it in the table of sharing requests.

The following parameters can be changed:

- The request's name
- The assays requested
- Institutions or locations can be added but not removed.

The following parameters cannot be changed:

- The start date
- Any of the data-aggregation options

### Note

After a sharing request has been configured to apply to institutions or to locations in the **Applies to** pane, that setting cannot be changed.

## 4.3 Delete a Sharing Request

A supranational institution can delete a sharing request it has created whether or not the receiving institutions have accepted the request. Follow these steps to delete a request.

1. Open the request by clicking it in the table of sharing requests.  
The sharing request's window opens.

### Note

Only people with the Supranational Institution Admin role can delete a sharing request.

2. Click **Delete** just above the right side of the **Properties** pane.
3. When the confirmation dialog opens click **Yes**.

The sharing request is deactivated and removed from the table of requests. Any data received from the sharing institution is removed from the requesting institution's aggregated data. (Deleting a sharing request does not remove data from the sending institution's database.) The request also is removed from the table of requests at the sharing institution. This occurs whether or not the institution has accepted the request.

## 4.4 Responding to Sharing Requests

Institutions that receive sharing requests can follow these steps to address the requests.

1. Under **ADMIN**, click the arrow next to the **Configurations** item to expand it. Then click **Sharing Requests**.

The list of pending and accepted sharing requests for your institution opens in a table. The table can be sorted by status, alphabetically by each request's name, or by the institution making the request.

Status	Name	Institution	Assays	Periodicity	Data category	Granularity	Data from
PENDING	BCR 2020s overview	CEP Supranational Medical East	Xpert BCR-ABL Monitor, Xpert BCR-ABL Ultra, Xpert BCR-ABL Ultra p190	Monthly	None	Laboratory	1/1/20
PENDING	C.Diff Metro Portland	CEP Supranational Medical East	Xpert C. difficile	Daily	None	Test Result with Patient Demographics	1/1/22
PENDING	Carba-R West	CEP Supranational Medical East	Xpert Carba-R	Real time	None	Laboratory	1/1/21
ACCEPTED	Multnomah schools	CEP Supranational Medical East	Xpert Norovirus	Daily	None	Test Result with Patient Demographics	1/1/21
ACCEPTED	Recent STD incidence	CEP Supranational Medical East	Xpert HPV 16,18-45, Xpert HPV HR, Xpert HIV Qual, Xpert HIV Quant	Daily	None	Subdistrict	1/1/22
PENDING	Recent Throat and Lung	CEP Supranational Medical East	Xpert Flu/RSV XC (Flu only), Xpert Flu, Xpert MRSA/SA Nasal Complete, Xpert MRSA, Xpert Xpress SARS-CoV-2/Flu/RSV, Xpert Xpress Strep A	Daily	None	Test Result with Patient Demographics	1/31/22

Items per page: 25    1 - 6 of 6    |< < > >|

2. Click a request in the table to open it. Review what details the requesting supranational institution has asked for.

Properties

**Requested by** CEP Supranational Medical East

**Requested data since** 1/1/21, 4:00 PM

**Involved assays** Xpert Norovirus

---

Requested data aggregation options

**Granularity** Test Result with Patient Demographics

**Time aggregation** Daily

**Delay** 1 days

---

Pending requests

**Assays** Xpert Norovirus

**Status** PENDING

**Share request data consent**

User **CEPSupra InstAdminOne** of Institution **CEP Supranational Medical East** has created a sharing request to view data from your organization. The granularity of this sharing request is **Test Result with Patient Demographics**, meaning **Test Result with Patient Demographics explanation**. Please accept or reject the request.


3. At the bottom of the sharing request screen is a box called **Share request data consent**. Click either **Agree to share this requested data** or **Reject sharing this requested data**.

**Important** The institution that receives a sharing request is responsible for complying with applicable privacy and data-protection laws.

## 4.5 Stopping Data Sharing

An institution that has accepted a sharing request can follow these steps to stop sharing data.

1. Open the request by clicking it in the table of sharing requests.  
The sharing request's window opens.
2. In the **Data being shared** pane at the bottom of the screen click **Stop sharing**.



Cepheid C360 stops sending new data to the requesting institution. Data previously shared with the requesting institution is no longer visible to the supranational institution. Stopped sharing requests remain in the list of requests.

## 4.6 Restarting Stopped Data Sharing

After having stopped a sharing request, you can follow these steps to start sharing data again.

1. Open the request by clicking it in the table of sharing requests.  
The sharing request's window opens. The **Stop sharing** button in the **Data being shared** pane now reads **Start sharing this data**.
2. Click **Start sharing this data**.  
Cepheid C360 begins sending the requested data to the supranational institution that initiated the request. If the request's start date was before or during the period the sharing request was stopped, any missed data is retrieved retroactively.

